Summary of Safety and Clinical Performance

FemoSeal™ Vascular Closure System

Manufacturer	Terumo Medical Corporation 265 Davidson Ave, Suite 320 Somerset, NJ 08873 USA
Manufacturer's SRN	US-MF-000019594
Management Representative/PRRC	Kimberly Feitl Vice President, Quality
Basic UDI-DI	38970FS79
EMDN number	C900199 – Haemostasis Systems-Other
Device Class	Class III
Year Device was First CE Marked	2005
EU Representative	Terumo Europe, N.V. Interleuvenlaan 40, 3001 Leuven Belgium SRN: BE-AR-000001433
Notified Body	NSAI CE 0050

Table of Contents

1.0	Summary for Patients	4
1.1	Intended Use of the Device.....	5
1.1.1	Intended Purpose.....	5
1.1.2	Indications for Use.....	5
1.1.3	Contraindications.....	5
1.1.4	Intended patient population.....	5
1.2	Device Description.....	6
1.2.1	Description of the Device.....	6
1.2.2	Raw Materials & Components.....	7
1.3	Risks and Warnings.....	8
1.3.1	Residual Risks and Undesirable Effects.....	9
1.3.2	Warnings and Precautions.....	10
1.3.3	Field Safety Corrective Actions.....	10
1.4	Summary of Clinical Evaluation and Post-Market Clinical Follow-Up.....	10
1.4.1	Clinical Background of the Device.....	10
1.4.2	Clinical Evidence Supporting CE Marking.....	11
1.4.3	Safety.....	23
1.5	Possible Therapeutic Alternatives.....	23
2.0	Bibliography	25
3.0	Revision History	28

This Summary of Safety and Clinical Performance (SSCP) is intended to provide public access to an updated summary of the main aspects of the safety and clinical performance of the FemoSeal™ Vascular Closure System (VCS).

A summary of the safety and clinical performance of the device, intended for patients, is provided in **Section 1.0**.

1.0 SUMMARY FOR PATIENTS

The information presented below is intended for patients or lay persons.

This SSCP is not intended to give general advice on the treatment of a medical condition. Please contact your healthcare professional in case you have questions about your medical condition or about the use of the device in your situation.

This SSCP is not intended to replace an Implant Card or the Instructions for Use to provide information on the safe use of the device.

Manufacturer	Terumo Medical Corporation 265 Davidson Ave, Suite 320 Somerset, NJ 08873 USA
Manufacturer's SRN	US-MF-000019594
Management Representative/PRRC	Kimberly Feitl Vice President, Quality
Basic UDI-DI	38970FS79
EMDN number	C900199 – Haemostasis Systems-Other
Device Class	Class III
Year Device was First CE Marked	2005
EU Representative	Terumo Europe, N.V. Interleuvenlaan 40, 3001 Leuven Belgium SRN: BE-AR-000001433
Notified Body	NSAI CE 0050

1.1 Intended Use of the Device

1.1.1 Intended Purpose

FemoSeal™ Vascular Closure System is a medical device used to close a hole in an artery (common femoral arterial puncture) of your groin.

1.1.2 Indications for Use

FemoSeal™ Vascular Closure System is approved for use in closing a hole in the artery in the groin.

1.1.3 Contraindications

FemoSeal™ is not to be used (contraindicated) in patients with holes in the artery that are larger than 2.33 mm (7 French) diameter.

1.1.4 Intended patient population

FemoSeal™ is approved to be used in patients 18 years of age and older.

The use of FemoSeal™ has not been tested in the following patients:

- Patients whose body's defense system attacks itself (auto-immune disease).
- Patients getting blood thinning medication (therapeutic thrombolysis).
- Patients with narrow arteries (clinically significant peripheral vascular disease) at the puncture site.
- Patients with high blood pressure (Uncontrolled hypertension of > 220 mmHg first number in a blood pressure reading [systolic] or > 110 mmHg second number in a blood pressure reading [diastolic]).
- Patients with a bleeding disorder (thrombocytopenia) or low iron (anemia).
- Patients with small arteries (common femoral artery smaller than 5mm).
- Patients with kidney or liver disease (renal and/or hepatic impairment).
- Patients with other diseases (relevant co-morbidity).
- Patients of a specific race or from a group that shares for example a common language or traditions (ethnicity).
- Patients who had a heart attack (myocardial infarction) within 72 hours.
- Patients with implanted device at the hole in the artery (Stent or vascular graft).
- Patients who are having a baby or breast feeding.


- Patients below the age of 18 years.

1.2 Device Description

1.2.1 Description of the Device

FemoSeal™ is a device used to close a hole in the femoral artery (arterial puncture). FemoSeal™ implantable parts include two degradable plastic (polymer) discs that are left in the body and a suture. The discs make a sandwich like seal held together by suture (multifilament). See **Figure 2.1** and **Figure 2.2** below.

Figure 2.1 FemoSeal™ VCS Unit with the Guidewire and Dilator

1. Guidewire 0.038" (0.97 mm)

2. Artery stretching tool (Dilator)


3. FemoSeal™ Unit

4. Tube that allows passage in and out of artery (Sheath)

Figure 2.2 Inner Seal and Outer Locking Disc, Respectively

Figure 2.3 Inner Seal and Outer Locking Disc Deployed

The discs left in the body dissolve within eighteen (18) months. The suture (multifilament) dissolves within two (2) to three (3) years.

1.2.2 Raw Materials & Components

FemoSeal™ parts that are left in the body will dissolve (degrade) over time. Pictures in **Figure 2.3**. **Table 2.1** below lists the materials in the device.

Table 2.1 Raw Materials in FemoSeal™ VCS

Part Number	Description	Materials
1	0.038" (0.97 mm) Guidewire with a guidewire J-straightener	Guidewire: 304 Stainless Steel J-Straightener: Polypropylene, purple pigment Clip: Polyethylene Tubing: Polyethylene
2	FemoSeal Dilator	Hub: Tetrahydrofuran, Blue pigment, Polybutylene terephthalate (PBT) Tube: High density polyethylene (HDPE), black pigment Lubricant: Silicone Dow Corning 360
3/4	FemoSeal Unit	Molded RD7 – Inner Seal and Outer locking Disc * Copolymer between glycolide, trimethylene carbonate, ε-caprolactone, and TMP: trimethylolpropane (initiator) Multifilament * Segmented copolymer between L,L-Lactide, trimethylene

			carbonate, ε-caprolactone, and 1,3 propanediol (initiator) Coating is copolymer between glycolide, ε-caprolactone, and L-lysine
		Tamping Tube	Polypropylene
		Pusher	304V Stainless steel
		Slider	Polybutylene terephthalate (PBT)
		Housing	Polybutylene terephthalate (PBT)
		Button	Polybutylene terephthalate (PBT)
		Button Lid	Polyoxymethylene (POM)/Polyacetal Copolymer
		Sleeve	Polybutylene terephthalate (PBT),
		Spacer	Polybutylene terephthalate (PBT)
		Cone	Polypropylene
		Housing Lid	Polybutylene terephthalate (PBT)
		Safety Catch	Polybutylene terephthalate (PBT), blue pigment
		Sleeve Lid	Polybutylene terephthalate (PBT)
		Tube Gasket	Silicone shore A 70
		Housing Gasket	Silicone shore A 50
		Spring	302 Stainless Steel
		Cone Housing Sheath	Polypropylene Silicone Dow Corning 360, Hexane

* Implantable portion of device; implantable components are MRI safe

1.3 Risks and Warnings

Contact your healthcare professional if you believe that you are experiencing side effects related to the device or its use or if you are concerned about risks. This document is not intended to replace a consultation with your healthcare professional if needed. Your doctor has considered the risks and decided that this device is appropriate for your procedure.

1.3.1 Residual Risks and Undesirable Effects

Harms associated with the use of the FemoSeal™ VCS as a result of residual risks are listed in **Table 2.2**. The occurrence rate is based on data collected from September 01, 2019 through August 31, 2024.

Table 2.2. FemoSeal Adverse Events/Complications

Known or Foreseeable Harm	Occurrence Rate (%)	Mitigating Factors
Allergic reaction	0.0000	Actions and plans to reduce risk as far as possible have been addressed by product design features, manufacturing guidelines, product labeling, and physician training.
Weakening of the artery wall (Aneurysm)	0.0000	
Abnormal connection between an artery and a vein (AV Fistula)	0.0001	
Blood loss / bleeding	0.0280	
Death	0.0002	
Bleeding under the skin (Ecchymosis)	0.0000	
Blockage of an artery (Embolism)	0.0000	
Irritation due to the device (Foreign body reaction)	0.0000	
Blood Loss (Hematoma)	0.0027	
Release of blood from a broken blood vessel (Hemorrhage)	0.0001	
Infection	0.0002	
Irritation, reddened area (Inflammation)	0.0000	
Loss of sensation (Numbness)	0.0000	
Pain	0.0000	
Patient discomfort	0.0004	
Procedure delay	0.0065	
An incomplete healing or sealing of the artery where the hole was	0.0004	

Known or Foreseeable Harm	Occurrence Rate (%)	Mitigating Factors
(Pseudoaneurysm/false aneurysm)		
Bleeding into the abdominal space (Retroperitoneal Bleed)	0.0002	
Infection (Sepsis)	0.0000	
Blockage of an artery due to a clot (Thromboembolism)	0.0000	
Blockage of an artery due to a clot (Thrombosis)	0.0001	
Blockage of an artery (Vessel occlusion or lower limb ischemia)	0.0002	
Tear in the artery (Vessel perforation)	0.0000	
Tear in the artery (Vessel tissue dissection/laceration)	0.0003	

1.3.2 Warnings and Precautions

- Follow doctor’s orders.
- Call your doctor if you have bleeding, pain, or signs of infection. Signs of infection are redness, fever, discharge and/or warmth at the site. Signs of a hematoma are swelling, lump in groin.
- Carry the Patient Information Card for the next 18 months.

1.3.3 Field Safety Corrective Actions

There was one recall involving FemoSeal™. The voluntary recall was due to an increase in product complaints. The failure was caused during manufacturing and has been corrected.

1.4 Summary of Clinical Evaluation and Post-Market Clinical Follow-Up

1.4.1 Clinical Background of the Device

FemoSeal™ has been on the market for 19 years.

1.4.2 Clinical Evidence Supporting CE Marking

FemoSeal™ approval is based on a study of 80 patients. The study showed that FemoSeal™ VCS was able to safely close the artery in at least 95% of patients.

Table 2.3 is a list of studies with FemoSeal™ VCS. The original version of the device was on the market until 2006. The current device version is in the marketplace today.

Table 2.3 Summary of Pivotal Studies

Pre-Market Functional Animal Testing	Device Version
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vasculature System. Functional Implantation Test in Sheep 12-, 15-, and 18-month Time Periods. (R1051-1), 2000, Radi Medical System AB (performed by Biomatech, France).	11200 (original device)
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep 12-, 24-, and 36-Week Time-Periods. (R1050-1), 2002, Radi Medical System AB (performed by Biomatech, France).	11200 (original device)
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep and in the Pig for 5 Weeks, (R1052-01) 2002, Radi Medical System AB (performed by Biomatech, France).	11200 (original device)
First In Human Clinical Study	
Multi-Centre Clinical Trial of the FemoSeal Vascular Closure System for Sealing Femoral Arterial Punctures After Diagnostic/Interventional Cardiology Procedures (EU-SE-001); 2004, Radi Medical Systems AB.	11200 (original device)
Completed Post-Market Clinical Studies	
PMS FemoSeal® VCS #11202 (R1756-01) (2006)	C11202 (current device)
Confirmatory Assessment of FemoSeal™ Vascular Closure System in 7F-Sheath. (2008)	C11202 (current device)
Post-Marketing Surveillance (PMS) Study FemoSeal® VCS #11202 (R1756-03) (2009)	C11202 (current device)
The FemoSeal Vascular Closure System (VCS) Registry: A Prospective, Multi-Center, Observational Study in Europe (2022)	C11202 (current device)
Planned / Ongoing Post-Market Clinical Studies	
FemoSeal™ Vascular Closure System PMCF Clinical Survey Protocol: Addressing Gaps Identified in CER-FS-2023	C11202 (current device)
Post Market Clinical Follow Up Activities	
20210292 - FemoSeal Vascular Closure System Clinical Survey Report (2020)	C11202 (current device)

20210298 – Guidewire Clinical Survey Report (2020)	C11202 (current device)
FemoSeal™ Vascular Closure System PMCF Clinical Survey Report: Addressing Gaps Identified in CER-FS-2021 and CER-FS-2022 (2024)	C11202 (current device)

Table 2.4 below describes all pre-market studies for FemoSeal™ VCS.

Table 2.4 Summaries of Pre-market Studies

Pre-Market Functional Animal Testing					
Name of Study	Year / CI	Study Type	Objectives	Milestones	Potential Acceptance Criteria
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vasculature System. Functional Implantation Test in Sheep 12-, 15-, and 18-month Time Periods. (R1051-1), 2000, Radi Medical System AB (performed by Biomatech, France).	2000 Radi Medical System AB Performed by Biomatech, France	Animal study; Exploratory investigation	The aim of the study was to evaluate the biocompatibility and biofunctionality of a biodegradable implant designed to achieve hemostasis in the vascular system. Report addresses reserve animals from Study Number 862 (Document # 1051-01)	N = 3 adult female sheep, 1 animal per time point Follow up for histopathologic analysis, macroscopic grading, and ultrastructural analysis <ul style="list-style-type: none"> 12 months 15 months 18 months Color doppler ultrasound follow up – 30 days at carotid and femoral arterial access sites	Arterial patency <ul style="list-style-type: none"> Post deployment blood flow assessment – 2/3 blood flow reduction associated with temporary arterial vasospasm that was attributed to procedure, both occurred at left common carotid artery Color doppler ultrasound assessment (4 weeks) – 3/3 (100%) without major blood flow reduction 30 days after implantation Macroscopic observations <ul style="list-style-type: none"> 12 months – patency confirmed macroscopically by absence of anatomically visible occlusion; no macroscopic local intolerance lesions (hemorrhage, necrosis, or neovascularization); encapsulation of internal disc appeared to be marked 15 months - macroscopic observation showed no signs of inflammation; no residues of inner or outer discs identified; encapsulation tissue located in area of outer disc appeared the same as 12-month observation 18 months – macroscopic observation showed no signs of local intolerance, outer disc

					<p>visible (3/3 cases), degradation of inner disc complete; 1/3 (Left common carotid artery) – depression of artery wall observed at implant site (could be related to device deployment)</p> <p>Conclusions</p> <ul style="list-style-type: none"> • Macroscopically, 12, 15 and 18 months after implantation, no necrotic, degenerative, or thrombotic signs were detected. • The degradation process of the inner and outer discs appeared to progress between 12 and 18 months and seemed complete at 18 months, with a whitish tissue covering the implanted area. • Histologically, after 18 months, polymer material was no longer visible but an endoluminal fibrous tissue was present. No necrotic, degenerative, or thrombotic lesions were noted.
<p>Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep 12-, 24-, and 36-Week Time-Periods. (R1050-1), 2002, Radi Medical System AB (performed by Biomatech, France).</p>	<p>2002 Radi Medical System AB</p> <p>Performed by Biomatech, France</p>	<p>Animal study; Exploratory investigation</p>	<p>The aim of the study was to evaluate the biocompatibility of a biodegradable implant designed to achieve hemostasis in the vascular system.</p>	<p>Follow up (after implantation)</p> <ul style="list-style-type: none"> • 12 weeks • 24 weeks • 36 weeks <p>N = 12 sheep, 3 animals per time point, 3 reserve animals sacrificed at 12, 15, and 18 months</p> <p>Implantation sites – right and left femoral arteries, right and left common carotid arteries</p> <p>Color doppler ultrasound analysis</p>	<p><u>Patency of carotid and femoral arteries after 4 weeks</u></p> <ul style="list-style-type: none"> • Major blood flow reduction 0/12 (0%) <p><u>No major blood flow reduction</u></p> <ul style="list-style-type: none"> • Severe blood flow impairment related to device insertion 0/12 (0%) • Blood flow reduction in left common carotid 2/12 (16.6%) due to temporary arterial vasospasm (attributed to surgery, not device) <p><u>Blood flow impairment by follow up</u></p> <p>12 weeks – 0/12 (0%) 24 weeks – 0/12 (0%) 36 weeks – 0/12 (0%)</p>

					<p>Blood pressure measurements (3 sheep – 36 weeks)</p> <ul style="list-style-type: none"> No major blood flow reduction detected <p>Macroscopic data for integration and degradation (Grade 0-4)</p> <p><u>Carotid sites</u></p> <p><i>Internal side</i></p> <p>-encapsulation</p> <p>12 week 1.3</p> <p>24 week 0.8</p> <p>36 week 0.7</p> <p>-material degradation</p> <p>12 week 2.2</p> <p>24 week 2.3</p> <p>36 week 0.7</p> <p><i>External side</i></p> <p>-encapsulation</p> <p>12 week 2.0</p> <p>24 week 1.0</p> <p>36 week 0.8</p> <p>-material degradation</p> <p>12 week 2.0</p> <p>24 week 2.7</p> <p>36 week 3.2</p> <p><u>Femoral sites</u></p> <p><i>External side</i></p> <p>-encapsulation</p> <p>12 week 1.0</p> <p>24 week 2.7</p> <p>36 week 1.0</p> <p>-material degradation</p> <p>12 week 2.0</p> <p>24 week 2.4</p> <p>36 week 3.0</p> <p><i>Internal side</i></p> <p>-encapsulation</p> <p>12 week 2.0</p> <p>24 week 1.2</p> <p>36 week 1.0</p> <p>-material degradation</p> <p>12 week 1.0</p> <p>24 week 2.2</p> <p>36 week 3.0</p>
--	--	--	--	--	--

					<p><u>Histological data for integration, degradation, and inflammatory parameters (mean of relevant values, Index 0-4) n=6</u></p> <p><u>Carotid sites</u></p> <p>External fibrosis 12 weeks – 2.3 24 weeks – 2.3 36 weeks – 1.8</p> <p>Internal neointima 12 weeks – 3.0 24 weeks – 2.7 36 weeks – 3.3</p> <p>Material degradation 12 weeks – 1.7 24 weeks – 3.0 36 weeks – 3 +</p> <p>Inflammatory parameters (macrophages) 12 weeks – 3.0 24 weeks – 1.7 36 weeks – 1.7</p> <p><u>Femoral sites</u></p> <p>External fibrosis 12 weeks – 1.3 24 weeks – 1.5 36 weeks – 1.7</p> <p>Internal neointima 12 weeks – 3.7 24 weeks – 3.3 36 weeks – 3.8</p> <p>Material degradation 12 weeks – 2.0 24 weeks – 3.0 36 weeks – 3+</p> <p>Inflammatory parameters (macrophages) 12 weeks – 3.5 24 weeks – 1.8 36 weeks – 1.8</p> <p><u>Histomorphological patency values</u></p> <p>Carotid sites 12 week 74.0% 24 week 81.9% 36 weeks -89.43%</p> <p>Femoral sites 12 week – 42.2%</p>

					24 week – 73.21% 36 weeks – 78.71% <u>US vessel patency (mean %)</u> Carotid site 12 week 94.6% 24 week – 96.7% 36 week – 94.7% Femoral sites 12 week – 66.9% 24 week – 63.6% 36 week – 75.1% Conclusions <ul style="list-style-type: none"> No significant local intolerance sign was detected (absence of visible inflammatory, necrotic, or degenerative lesions) at 12, 24, or 36 weeks No significant adverse tissular reaction recorded with regard to host and device Degradation of implant was nearly complete at 36 weeks
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep and in the Pig for 5 Weeks, (R1052-01) 2002, Radi Medical System AB (performed by Biomatech, France).	2002 Radi Medical System AB Performed by Biomatech, France	Animal study	The purpose of the study was to evaluate the biocompatibility and biofunctionality of a biodegradable implant designed to achieve haemostasis in the vascular system. The tested device was intended to enter the puncture hole in arteries to achieve haemostasis. One disc entered the puncture hole in the artery and the other tightened the internal disc from outside the artery. Both discs were inserted through an introducer sheath. The discs were designed to prevent bleeding,	Implantation period / follow up is 5 weeks for all animals 4 sheep, 1 pig 14 closure elements 8 closure elements implanted in 2 sheep, 2 closure elements implanted in 1 pig When two closure elements were implanted into same femoral artery, they were implanted at a distance > 30 mm	Any adverse signs criteria (inflammation, necrosis, hemorrhage, or any other lesion) or any thin pseudo-intimal layer formation were recorded. Particular attention was devoted to initial tissular integration of the device, with a special emphasis on the internal disc encapsulation. Macroscopically, no significant signs of local intolerance or thrombus formation were observed. The macroscopic observation did not reveal any difference between the sheep or the pig after implantation of the closure elements in the femoral arteries. Macroscopically,

			<p>to heal into the vessel wall and finally degrade and be absorbed. The study was designed to collect and analyze data</p> <p>after functional implantation in the sheep and pig to evaluate the capacity of the implant to maintain physiologic function</p> <p>when used in the circulatory system, to determine the response of the host and the response of the device (physical integrity, tissular encapsulation at sacrifice) and to evaluate the degradation of the implant</p>		<p>no significant signs of local intolerance or thrombus formation were observed. The macroscopic observation did not reveal any difference between the sheep or the pig after implantation of the closure elements in the femoral arteries.</p> <p>The histopatological analysis of the implantation site in sheep showed no local adverse reactions. The inner and outer discs showed signs of degradation and were in some cases broken into smaller parts. The inner discs were completely integrated into neointimal tissue. In the pig the histopatological analysis showed slightly more inflammatory signs around the inner and outer discs. In both animals, the multifilament exhibited slightly more inflammatory signs than in the disc material, which is normal for braided sutures due to the larger surface area.</p> <p>The vessel patency was measured after 5 weeks using three different methods, histomorphometric, histomorphometric in combination with measurement of vessel diameter by use of caliper and by the use of color doppler ultrasonography. The mean vessel patency in sheep was between 60</p>
--	--	--	--	--	--

					and 80 % depending on method (12 observations). In the pig, the mean vessel patency was between 40 and 80 % (2 observations).
First In Human Clinical Study					
Name of Study	Year / CI	Study Type	Objectives	Milestones	Potential Acceptance Criteria
Multi-Centre Clinical Trial of the FemoSeal Vascular Closure System for Sealing Femoral Arterial Punctures After Diagnostic/Interventional Cardiology Procedures (EU-SE-001); 2004, Radi Medical Systems AB.	2004 Radi Medical Systems AB	Open label, single arm, multi-centre trial	The objectives of the study were to determine the safety, performance, and effectiveness of the FemoSeal Vascular Closure System (FemoSeal System) in general, and specifically the Closure Elements, in patients undergoing cardiac diagnostic or interventional cardiology procedures using femoral artery access.	2-part study Total n = 80 Part 1 n = 20 Part 2 n = 60 Indications: Part 1 – diagnostic cardiac catheterization using femoral artery as site of access Part 2 – patients undergoing diagnostic and interventional cardiac procedures using femoral artery as the site of access 30 day follow up	Outcomes <u>hemostasis achieved with FS without compression</u> - Part 1 20/20 (100%) - Part2 58/60 (97%) Median time to hemostasis at time of procedure (part 2) – 1.0 minutes (Range immediate – 8 minutes) <u>Time to ambulation (minutes)</u> Part 1 – Median: 33.0 mean: 67.8 Part 2 – median: 45 mean: 87.5 <u>Incidence of major vascular complications – 30 day follow up</u> Part 1 -Free of vascular complications 20/20 (100%) Part 2 Pseudoaneurysm - 1/59 (1.69%) <u>Treatment failure</u> Part 1 – 0/20 (0%) Part 2 – 3/60 (5%)

					<ul style="list-style-type: none"> 2 due to the failure to achieve hemostasis without use of compression 1 due to major vascular complication, pseudoaneurysm <u>Persistent oozing</u> Part 1 – 0/20 (0%) Part 2- 8/60 (13%) <u>Hematoma</u> Part 1 – 0/20 (0%) Part 2 – 6/60 (10%)
--	--	--	--	--	---

Table 2.5 below summarizes all studies completed once the device went to market.

Table 2.5 Post-Market Studies

Completed Post-Market Clinical Studies					
Name of Study	Year / CI	Study Type	Objectives	Milestones	Potential Acceptance Criteria
PMS FemoSeal® VCS #11202 (R1756-01) (2006) User Acceptability of New Packaging (SIS) (2006)	2006 RADI	Post Market Study, User Acceptability of New Packaging	Post-market study carried out on the model update FemoSeal Vascular® Closure System #11202 which has a modified delivery system compared to the first model, #11200, to document functionality. The objective of the study was to get feedback on the revised Instructions for Use (IFU), the training project for physicians and nurse, to confirm the assessments made in the risk analysis and to give input to any product quality improvements.	100 patients <ul style="list-style-type: none"> 63 male / 37 females Mean age 68 years for female, 65 years for males 2 Swedish hospitals, 6 physicians 112 devices Data captured on patient record forms (CRF) Planned observation time – 30 minutes from time of application	Outcomes 1. Illustrations/Instructions in IFU considered informative and adequate 2. Immediate hemostasis (91/100)
Confirmatory Assessment of FemoSeal™ Vascular Closure System in 7F-Sheath. (2008)	Dr. Nicolas Moes Innsbruck Medical	Post -Market Confirmatory Assessment	The objective of this study was to assess safety and performance of FemoSeal® used during normal clinical routine in	50 consecutive patients Indication – diagnostic or invasive cardiac catheterization from common femoral artery	Primary endpoints <ul style="list-style-type: none"> Primary efficacy endpoint – TTH (time to hemostasis) with or without oozing but without any additional

	University, Austria		<p>patients undergoing diagnostic cardiac catheterization or percutaneous cardiac interventional procedure performed with 7F-sheath.</p>	<table border="1" style="font-size: small;"> <tr><td>mean age</td><td>62.8 ± 11.1</td></tr> <tr><td>men</td><td>35 (70%)</td></tr> <tr><td>women</td><td>15 (30%)</td></tr> <tr><td>PCI cases</td><td>15</td></tr> <tr><td>age (yr)</td><td>62.8 ± 11.1</td></tr> <tr><td>body mass index (kg/m²)</td><td>25.4 ± 6.7</td></tr> <tr><td>systolic BP (mmHg)</td><td>131.2 ± 20.9</td></tr> <tr><td>diastolic BP (mmHg)</td><td>77.9 ± 10.6</td></tr> <tr><td>diabetes mellitus</td><td>9 (18%)</td></tr> <tr><td>arterial hypertension</td><td>34 (68%)</td></tr> <tr><td>previous coronary intervention</td><td>13 (26%)</td></tr> <tr><td>preprocedure aspirin</td><td>50 (100%)</td></tr> <tr><td>preprocedure clopidogrel</td><td>50 (100%)</td></tr> <tr><td>Procedural unfractionated heparin</td><td>15 (30%)</td></tr> <tr><td>Procedural abciximab</td><td>2 (4%)</td></tr> </table>	mean age	62.8 ± 11.1	men	35 (70%)	women	15 (30%)	PCI cases	15	age (yr)	62.8 ± 11.1	body mass index (kg/m ²)	25.4 ± 6.7	systolic BP (mmHg)	131.2 ± 20.9	diastolic BP (mmHg)	77.9 ± 10.6	diabetes mellitus	9 (18%)	arterial hypertension	34 (68%)	previous coronary intervention	13 (26%)	preprocedure aspirin	50 (100%)	preprocedure clopidogrel	50 (100%)	Procedural unfractionated heparin	15 (30%)	Procedural abciximab	2 (4%)	<p>compression, manual or mechanical.</p> <p>Safety endpoint - incidence of adverse events - minor and major vascular complications and all other adverse events – serious and non-serious during duration of the study (i.e., until patient discharge).</p> <p>Outcomes</p> <ul style="list-style-type: none"> • Mean time to hemostasis • Diagnostic - 56.6 ± 29 seconds • Interventional 59.9 ± 16.2 seconds • All patients 57.8 ± 26.3 • All patients could be ambulated after 4 hours bed rest without any problems. <p>Safety outcomes</p> <ul style="list-style-type: none"> • 49/50 (98%) free of vascular complications until hospital discharge • Mean time to discharge 2.6 days • 1/50 (2%) minor bleeding next day after uneventful ambulation on day of angiography
mean age	62.8 ± 11.1																																		
men	35 (70%)																																		
women	15 (30%)																																		
PCI cases	15																																		
age (yr)	62.8 ± 11.1																																		
body mass index (kg/m ²)	25.4 ± 6.7																																		
systolic BP (mmHg)	131.2 ± 20.9																																		
diastolic BP (mmHg)	77.9 ± 10.6																																		
diabetes mellitus	9 (18%)																																		
arterial hypertension	34 (68%)																																		
previous coronary intervention	13 (26%)																																		
preprocedure aspirin	50 (100%)																																		
preprocedure clopidogrel	50 (100%)																																		
Procedural unfractionated heparin	15 (30%)																																		
Procedural abciximab	2 (4%)																																		
<p>Post-Marketing Surveillance (PMS) Study FemoSeal® VCS #11202 (R1756-03) (2009)</p>	2009 RADI	Post Market Study	<p>The study should inform if the customers/users perceived any problems with the new package, especially the non-sterile outer side of the Tuvek/PET-PE pouch.</p> <p>The study should also give information about if the treating physician perceives any differences in the application/deployment steps of the device compared to previously used devices.</p>	<p>25 devices</p> <p>4 Swedish hospitals</p> <p>16/25 (64%) of patients were PCI patients with heavy anti-coagulation medications</p>	<p>Outcomes</p> <ol style="list-style-type: none"> 1. Package related problems – 21/25 (84%) applications were successful; 4/25 (16%) were unsuccessful (1 - testing new introducer, 1 inner disc when through arteriotomy at backing step, button came up and outer disc was placed in incision tissue canal, 2 immediate hemostasis was not achieved despite no deployment issues.) 2. Immediate hemostasis 19/25 (76%) patients 3. Bleeding 11/25 (44%) – 9/25 (36%) oozing/bleeding, 2/25 (8%) hematomas <p>Total event rate 19%</p>																														

			The study should also give input to any further needs for product quality improvement.	
<p>The FemoSeal Vascular Closure System (VCS) Registry: A Prospective, Multi-Center, Observational Study in Europe (T138E4) (2022)</p>	<p>2022/Terumo Medical Corporation</p>	<p>Prospective, Multi-center, Observational, Post-market</p>	<p>To further demonstrate the safety and effectiveness of the FemoSeal™ VCS in achieving hemostasis of common femoral artery (CFA) access site in real-world subjects undergoing percutaneous endovascular procedures</p>	<ul style="list-style-type: none"> • Patient follow-up at 30 days (\pm 7 days) by hospital visit or telephone call. • Enrollment: December 2021 through July 2022. • 230 enrolled subjects undergoing diagnostic or interventional endovascular procedures in which FemoSeal™ is used. <p>3 European sites: France (92 patients), Germany (38 patients), Belgium (100 patients).</p>
				<ul style="list-style-type: none"> • Demographics/comorbidities: <ul style="list-style-type: none"> - Gender: 161 males (70 %)/69 females 69 (30%). - Mean age: 70 \pm 12. - BMI: 26.1 \pm 4.95 kg/m². - Hypertension: 160 patients (69.6 %). - Dyslipidemia: 128 (55.7 %). - Coronary artery disease: 42 (18.3%). - Current smoker: 86 (37.6 %). - Past smoker: 66 (28.8%). - Coronary heart disease: 34 (14%). - Previous myocardial infarction: 18 (7.8%). - Peripheral artery revascularization: 12 (5.2 %) carotid and 83 (37.5%) and low limb artery. - Antiplatelet therapy (at baseline): 211 (91,7 %), including 21 (9.1%) dual antiplatelet therapy (DAPT). - Oral anticoagulation: 40 (17.4%) patients. • Intervention characteristics: <ul style="list-style-type: none"> - Setting: inpatient clinic 162 cases (70.4%), outpatient clinic 68 cases (29.6%). - Prevailing indication: claudication 147 cases (63.9%). - Approach: 35 antegrade (15.3%), 194 retrograde (84.7%). • Primary Endpoints: <ul style="list-style-type: none"> - Combined safety (freedom from major complications of the access site limb within 6 hours post-

					<p>procedure) and effectiveness (successful puncture site hemostasis) endpoint achieved in 215/226 (95.1%) [95% CI: 91.46; 97.55] patients in the Full Analysis Set and 215/230 (93.5%) [95% CI 89.47, 96.30] in the FemoSeal™ Treated Set.</p> <ul style="list-style-type: none"> - Effectiveness endpoint was achieved in 219/226 (96.9%) CI [93.7, 98.7] of patients. - Safety endpoint achieved in 220/230 (95.2%) [95% CI: 92.15, 97.90] patients. <p>• Secondary endpoints:</p> <ul style="list-style-type: none"> - Freedom from minor access site complications for six hours post-procedure achieved in 225/230 (97.8%) [95% CI: 95.00; 99.29]. - Freedom from major and minor access site complications from six hours to 30 days post-procedure achieved in 219/230 (95.2%) [95% CI:91.60; 97.59] of patients. - Median Time to hemostasis with FemoSeal™ VCD: 0.42 (0.25; 0.50) minutes. - Median Time to ambulation: 5.00 (4.54; 5.50) minutes. - Median length of stay in hospital: 23.98 (22.72, 25.00) hours. - FemoSeal™ VCS usability: ‘easy’ or ‘very easy’ to deploy in 98.7% of cases (‘easy’ to deploy in 21 (9.2%), ‘very easy’ to deploy in 205 (89.5 %)) with no or low resistance at deployment in 226 (98.3%) of cases.
--	--	--	--	--	---

					152/230 (66.1%) of patients treated using guidewire included with the FemoSeal™ kit (GW 0.038") and recommended in the FemoSeal™ without any reported device deficiencies or adverse events.
--	--	--	--	--	--

1.4.3 Safety

FemoSeal™ has a shorter time to stop the bleeding at the procedure site compared to manual pressure. FemoSeal™ and manual pressure have similar risks. The clinical benefits of using FemoSeal™ outweigh the risks. Your doctor has considered the risks and decided that this device is appropriate for your procedure. See Section 2.1.4 for patient population information. If you have an existing health condition, speak with your doctor.

There were 4 post-market clinical studies for FemoSeal™. The 4 studies included a 2006 study with 100 patients, a 2008 study with 50 patients, 2009 study with 25 patients, and a 2022 prospective study with 230 patients. These studies confirmed the safety and effectiveness of FemoSeal™. These studies also confirmed that FemoSeal™ performs as it is intended.

In addition, a survey of physicians was performed to gather information on the guidewire used with the device unit. The survey gathered information from physicians on what type of guidewire they use during the procedure.

Terumo Medical Corporation continues to collect information on the device after its released in market. This is to demonstrate long-term safety and performance of the device.

1.5 Possible Therapeutic Alternatives

Manual pressure is the standard method to stop bleeding. It is the use of external hand pressure, by a member of the hospital staff, to the hole in the artery. While this method works, this can be uncomfortable to patients. They often must lay flat for long periods of time, possibly leading to longer hospital stays. On average it takes 21 minutes for bleeding to stop during manual compression.

There are other ways to stop the bleeding, such as mechanical compression devices which assist manual compression. These devices press on the wound mechanically instead of staff members applying pressure. They do not speed up the time to stop the bleeding or help patients walk sooner but they do reduce the need for staff to apply pressure manually. The risk of problems from these devices are about the same as manual compression.

Another option is vascular closure devices like FemoSeal™ VCS. These devices help close the artery opening which can stop bleeding faster and allow patients to move around sooner. A recent study looked at 13 different types of closure devices and found on average they stop bleeding in about 5 minutes.

No matter which method is used some risks can occur. The most common include bleeding, small artery bulges (pseudoaneurysms), abnormal connections between blood vessels (arteriovenous fistulas) and infection. A large study involving 14,000 patients found that these risks were similar for all methods, although bleeding was slightly more common with vascular closure devices.

When considering other treatments, contact your doctor who can consider your medical history.

2.0 BIBLIOGRAPHY

1. Alabduljabbar M, Alhamdan B, Alabdrabalnabi M. Access site complications of manual compression versus closure devices after lower limb revascularization: a systematic review. *International Journal of Medicine in Developing Countries*. 2022;6(12):1572-1580. doi:10.24911/ijmdc.51-1666716459
2. Chugh Y, Bavishi C, Mojadidi MK, et al. Safety of transradial access compared to transfemoral access with hemostatic devices (vessel plugs and suture devices) after percutaneous coronary interventions: A systematic review and meta-analysis. *Catheter Cardiovasc Interv*. Aug 2020;96(2):285-295. doi:10.1002/ccd.29061
3. Ciprian Cacuci A, Krankenberg H, Ingwersen M, et al. Access Site Complications of Peripheral Endovascular Procedures: A Large, Prospective Registry on Predictors and Consequences. *J Endovasc Ther*. Oct 2021;28(5):746-754. doi:10.1177/15266028211025044
4. Faraj J, Daneshmand A, Mwiapatayi BP. Large spontaneous median sacral artery pseudoaneurysm masquerading as hip pain. Article. *J Vasc Surg Cases Innov Tech*. Jun 2020;6(2):224-227. doi:10.1016/j.jvscit.2020.02.012
5. Gabrielli R, Siani A, Smedile G, Rizzo AR, Accrocca F, Bartoli S. Carotid Artery Stenting versus Carotid Endarterectomy in Terms of Neuroprotection DW-MRI Detected and Neuropsychological Assessment Impairment. Article. *Ann Vasc Surg*. Jan 2024;98:68-74. doi:10.1016/j.avsg.2023.05.046
6. Goueffic Y, Picquet J, Schneider F, et al. A Randomized Trial Comparing Polymer Versus Suture- Based Vascular Closure Devices for Arterial Closure Following Lower-Limb Arterial Endovascular Revascularization. *Cardiovasc Intervent Radiol*. Dec 2021;44(12):1883-1892. doi:10.1007/s00270-021-02940-z
7. Ha TH, Lee SH, Park SJ, Chung HH. Could real-time ultrasonography guidance be useful for the effective deployment of FemoSeal in common femoral arteriotomy? *Ultrasonography*. Jul 2021;40(3):449-454. doi:10.14366/usg.20127
8. Hassu J, Meyborg M, Gebauer K, Makowski L, Malyar N. 51. Embolization of closure device particle as a rare case of new-onset claudication - Normal diagnostic findings despite limiting symptoms. A Diagnostic Challenge. Conference Abstract. *Vasa*. Nov 2022;51(Suppl 109):1-21. doi:10.1024/0301-1526/a001037
9. Kaszczewski P, Kozubek H, Ostrowski T, et al. Rapid Formation and Hybrid Treatment of a Large Superior Mesenteric Artery Aneurysm. Article. *Am J Case Rep*. Oct 29 2023;24:e939558. doi:10.12659/AJCR.939558

10. Kaszczewski P, Leszczynski J, Elwertowski M, Maciag R, Chudzinski W, Galazka Z. Combined Treatment of Multiple Splanchnic Artery Aneurysms Secondary to Median Arcuate Ligament Syndrome: A Case Study and Review of the Literature. Article. Am J Case Rep. Aug 11 2020;21:e926074. doi:10.12659/AJCR.926074
11. Kennedy SA, Rajan DK, Bassett P, Tan KT, Jaber A, Mafeld S. Complication rates associated with antegrade use of vascular closure devices: a systematic review and pooled analysis. J Vasc Surg. Feb 2021;73(2):722-730 e1. doi:10.1016/j.jvs.2020.08.133
12. Langouet Q, Martinez R, Saint-Etienne C, et al. Incidence, predictors, impact, and treatment of vascular complications after transcatheter aortic valve implantation in a modern prospective cohort under real conditions. Article. J Vasc Surg. Dec 2020;72(6):2120-2129 e2. doi:10.1016/j.jvs.2020.03.035
13. Mach M, Maciejewski K, Ostrowski T, Maciag R, Sajdek M, Galazka Z. Endovascular Treatment of a Bilateral, Ruptured Angiomyolipoma in a Patient With Tuberous Sclerosis Complex. Cureus. Aug 2024;16(8):e66200. doi:10.7759/cureus.66200
14. Mayer K, Gewalt S, Morath T, et al. Comparison of Vascular Closure Devices vs Manual Compression After Femoral Artery Puncture in Patients on Oral Anticoagulation - Post Hoc Analysis of the ISAR-CLOSURE Trial. The Journal of invasive cardiology. Sep 2021 2021;33(9):E709-E715.
15. Messana G, Ambrosi L, Moramarco LP, Cionfoli N, Maestri M, Quaretti P. Testicular artery originating from the inferior mesenteric artery: an alert for interventionalists – A case report. Article. Radiology Case Reports. 2021;16(9):2710-2713. doi:10.1016/j.radcr.2021.06.059
16. Moramarco LP, Capodaglio CA, Quaretti P, et al. Multivessel endovascular therapy for undiagnosed vascular type Ehlers-Danlos syndrome. Successful percutaneous transcatheter coil embolization of hepatic artery pseudoaneurysm with stenting of right renal and iliac arteries in emergency setting. BJR Case Rep. Dec 1 2020;6(4):20200025. doi:10.1259/bjrcr.20200025
17. Morosetti D, Chiochi M, Argiro R, et al. Endovascular treatment of calcific lesions of the common femoral artery using atherectomy device associated with scoring balloon angioplasty in diabetic patients with high "major amputation" risk. Vascular. Jun 2022;30(3):463-473. doi:10.1177/17085381211019244
18. Noory E, Bohme T, Krause L, et al. Evaluation of the MYNX CONTROL Arterial Closure System for Achieving Primary Hemostasis after Arterial Femoral Access Following Peripheral Arterial Interventions, Compared to the FemoSeal(TM) Closure System. J Clin Med. Aug 12 2023;12(16):5255. doi:10.3390/jcm12165255

19. Oliny A, Bosse C, Charbonneau P, Haulon S, Fabre D. Closure Device Migration: An Unusual Cause of Acute Limb Ischaemia Following a Simple Endovascular Procedure. Article. EJVES Vasc Forum. 2024;61:113-115. doi:10.1016/j.ejvsf.2024.05.006
20. Pane F, Coppola M, Borzelli A, et al. A singular case of massive urethrorrhagia solved by transarterial embolization. Article. Egyptian Journal of Radiology and Nuclear Medicine. 2020;51(1):168. doi:10.1186/s43055-020-00299-8
21. Pang N, Gao J, Zhang B, et al. Vascular Closure Devices versus Manual Compression in Cardiac Interventional Procedures: Systematic Review and Meta-Analysis. Cardiovasc Ther. 2022 2022;2022:8569188. doi:10.1155/2022/8569188
22. Tagliaferro FB, Orgera G, Mascagni L, et al. FemoSeal((R)) vascular closure device for antegrade common femoral artery access: Safety and technical notes. J Vasc Access. Jan 2020;21(1):79- 85. doi:10.1177/1129729819854593
23. Tan R, Sieunarine K. Superficial Femoral Artery Pseudoaneurysm as a Delayed Complication of TurboHawk Atherectomy. Annals of Vascular Surgery. Oct 2020 2020;68:568.e17-568.e21. doi:10.1016/j.avsg.2020.04.006

3.0 REVISION HISTORY

SSCP revision number	Date issued	Change description	Revision validated by the Notified Body
SSCP-FS-2022	January 20, 2023	Initial Issue, upon MDR CE certification	<input checked="" type="checkbox"/> Yes Validation language: English <input type="checkbox"/> No (only applicable for class IIa or some IIb implantable devices (MDR, Article 52 (4) 2nd paragraph) for which the SSCP is not yet validated by the NB) <input type="checkbox"/> Other (explain):
SSCP-FS-2023 – Revision 1	March 15, 2024	Updated for the period between January 1, 2019 through December 31, 2023.	<input checked="" type="checkbox"/> Yes Validation language: English <input type="checkbox"/> No (only applicable for class IIa or some IIb implantable devices (MDR, Article 52 (4) 2nd paragraph) for which the SSCP is not yet validated by the NB) <input type="checkbox"/> Other (explain):
SSCP-FS-2024 Revision 1		Updated for the period between September 1, 2019, through August 31, 2024.	<input type="checkbox"/> Yes Validation language: English <input type="checkbox"/> No (only applicable for class IIa or some IIb implantable devices (MDR, Article 52 (4) 2nd paragraph) for which the SSCP is not yet validated by the NB) <input type="checkbox"/> Other (explain):
SSCP-FS-2024 Revision 2	June 20, 2025	Revised based on NSAI queries for validation	<input checked="" type="checkbox"/> Yes Validation language: English

			<p><input type="checkbox"/> No (only applicable for class IIa or some IIb implantable devices (MDR, Article 52 (4) 2nd paragraph) for which the SSCP is not yet validated by the NB)</p> <p><input type="checkbox"/> Other (explain):</p>
--	--	--	---

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Функционално тестване върху животни преди пускане на пазара	Версия на изделието
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vasculature System. Functional Implantation Test in Sheep 12-, 15-, and 18-month Time Periods. (R1051-1), 2000, Radi Medical System AB (performed by Biomatech, France).	11200 (оригиналното изделие)
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep 12-, 24-, and 36-Week Time-Periods. (R1050-1), 2002, Radi Medical System AB (performed by Biomatech, France).	11200 (оригиналното изделие)
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep and in the Pig for 5 Weeks, (R1052-01) 2002, Radi Medical System AB (performed by Biomatech, France).	11200 (оригиналното изделие)
Първо клинично проучване при хора	
Multi-Centre Clinical Trial of the FemoSeal Vascular Closure System for Sealing Femoral Arterial Punctures After Diagnostic/Interventional Cardiology Procedures (EU-SE-001); 2004, Radi Medical Systems AB.	11200 (оригиналното изделие)
Завършени клинични проучвания след пускане на пазара	
PMS FemoSeal® VCS #11202 (R1756-01) (2006)	C11202 (настоящото изделие)
Confirmatory Assessment of FemoSeal™ Vascular Closure System in 7F-Sheath. (2008)	C11202 (настоящото изделие)
Post-Marketing Surveillance (PMS) Study FemoSeal® VCS #11202 (R1756-03) (2009)	C11202 (настоящото изделие)
The FemoSeal Vascular Closure System (VCS) Registry: A Prospective, Multi-Center, Observational Study in Europe (2022)	C11202 (настоящото изделие)

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2.0 LITERATURA

1. Alabduljabbar M, Alhamdan B, Alabdrabalnabi M. Access site complications of manual compression versus closure devices after lower limb revascularization: a systematic review. *International Journal of Medicine in Developing Countries*. 2022;6(12):1572-1580. doi:10.24911/ijmdc.51-1666716459
2. Chugh Y, Bavishi C, Mojadidi MK, et al. Safety of transradial access compared to transfemoral access with hemostatic devices (vessel plugs and suture devices) after percutaneous coronary interventions: A systematic review and meta-analysis. *Catheter Cardiovasc Interv*. Aug 2020;96(2):285-295. doi:10.1002/ccd.29061
3. Ciprian Cacuci A, Krankenberg H, Ingwersen M, et al. Access Site Complications of Peripheral Endovascular Procedures: A Large, Prospective Registry on Predictors and Consequences. *J Endovasc Ther*. Oct 2021;28(5):746-754. doi:10.1177/15266028211025044
4. Faraj J, Daneshmand A, Mwiapatayi BP. Large spontaneous median sacral artery pseudoaneurysm masquerading as hip pain. *Article. J Vasc Surg Cases Innov Tech*. Jun 2020;6(2):224-227. doi:10.1016/j.jvscit.2020.02.012
5. Gabrielli R, Siani A, Smedile G, Rizzo AR, Accrocca F, Bartoli S. Carotid Artery Stenting versus Carotid Endarterectomy in Terms of Neuroprotection DW-MRI Detected and Neuropsychological Assessment Impairment. *Article. Ann Vasc Surg*. Jan 2024;98:68-74. doi:10.1016/j.avsg.2023.05.046
6. Goueffic Y, Picquet J, Schneider F, et al. A Randomized Trial Comparing Polymer Versus Suture- Based Vascular Closure Devices for Arterial Closure Following Lower-Limb Arterial Endovascular Revascularization. *Cardiovasc Intervent Radiol*. Dec 2021;44(12):1883-1892. doi:10.1007/s00270-021-02940-z
7. Ha TH, Lee SH, Park SJ, Chung HH. Could real-time ultrasonography guidance be useful for the effective deployment of FemoSeal in common femoral arteriotomy? *Ultrasonography*. Jul 2021;40(3):449-454. doi:10.14366/usg.20127
8. Hassu J, Meyborg M, Gebauer K, Makowski L, Malyar N. 51. Embolization of closure device particle as a rare case of new-onset claudication - Normal diagnostic findings despite limiting symptoms. *A Diagnostic Challenge. Conference Abstract. Vasa*. Nov 2022;51(Suppl 109):1-21. doi:10.1024/0301-1526/a001037
9. Kaszczewski P, Kozubek H, Ostrowski T, et al. Rapid Formation and Hybrid Treatment of a Large Superior Mesenteric Artery Aneurysm. *Article. Am J Case Rep*. Oct 29 2023;24:e939558. doi:10.12659/AJCR.939558

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Funktionstest på dyr før markedsføring	Version af udstyret
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vasculature System. Functional Implantation Test in Sheep 12-, 15-, and 18-month Time Periods. (R1051-1), 2000, Radi Medical System AB (performed by Biomatech, France).	11200 (originalt udstyr)
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep 12-, 24-, and 36-Week Time-Periods. (R1050-1), 2002, Radi Medical System AB (performed by Biomatech, France).	11200 (originalt udstyr)
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep and in the Pig for 5 Weeks, (R1052-01) 2002, Radi Medical System AB (performed by Biomatech, France).	11200 (originalt udstyr)
Første kliniske forsøg med mennesker	
Multi-Centre Clinical Trial of the FemoSeal Vascular Closure System for Sealing Femoral Arterial Punctures After Diagnostic/Interventional Cardiology Procedures (EU-SE-001); 2004, Radi Medical Systems AB.	11200 (originalt udstyr)
Gennemførte kliniske forsøg efter markedsføring	
PMS FemoSeal® VCS #11202 (R1756-01) (2006)	C11202 (aktuel enhed)
Confirmatory Assessment of FemoSeal™ Vascular Closure System in 7F-Sheath. (2008)	C11202 (aktuel enhed)
Post-Marketing Surveillance (PMS) Study FemoSeal® VCS #11202 (R1756-03) (2009)	C11202 (aktuel enhed)
The FemoSeal Vascular Closure System (VCS) Registry: A Prospective, Multi-Center, Observational Study in Europe (2022)	C11202 (aktuel enhed)
Planlagte/igangværende kliniske forsøg efter markedsføring	
FemoSeal™ Vascular Closure System PMCF Clinical Survey Protocol: Addressing Gaps Identified in CER-FS-2023	C11202 (aktuel enhed)
Kliniske opfølgingsaktiviteter efter markedsføring	
20210292 - FemoSeal Vascular Closure System Clinical Survey Report (2020)	C11202 (aktuel enhed)
20210298 – Guidewire Clinical Survey Report (2020)	C11202 (aktuel enhed)
FemoSeal™ Vascular Closure System PMCF Clinical Survey Report: Addressing Gaps Identified in CER-FS-2021 and CER-FS-2022 (2024)	C11202 (aktuel enhed)

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2.0 BIBLIOGRAFI

1. Alabduljabbar M, Alhamdan B, Alabdrabalnabi M. Access site complications of manual compression versus closure devices after lower limb revascularization: a systematic review. *International Journal of Medicine in Developing Countries*. 2022;6(12):1572-1580. doi:10.24911/ijmdc.51-1666716459
2. Chugh Y, Bavishi C, Mojadidi MK, et al. Safety of transradial access compared to transfemoral access with hemostatic devices (vessel plugs and suture devices) after percutaneous coronary interventions: A systematic review and meta-analysis. *Catheter Cardiovasc Interv*. Aug 2020;96(2):285-295. doi:10.1002/ccd.29061
3. Ciprian Cacuci A, Krankenberg H, Ingwersen M, et al. Access Site Complications of Peripheral Endovascular Procedures: A Large, Prospective Registry on Predictors and Consequences. *J Endovasc Ther*. Oct 2021;28(5):746-754. doi:10.1177/15266028211025044
4. Faraj J, Daneshmand A, Mwiratayi BP. Large spontaneous median sacral artery pseudoaneurysm masquerading as hip pain. *Article. J Vasc Surg Cases Innov Tech*. Jun 2020;6(2):224-227. doi:10.1016/j.jvscit.2020.02.012
5. Gabrielli R, Siani A, Smedile G, Rizzo AR, Accrocca F, Bartoli S. Carotid Artery Stenting versus Carotid Endarterectomy in Terms of Neuroprotection DW-MRI Detected and Neuropsychological Assessment Impairment. *Article. Ann Vasc Surg*. Jan 2024;98:68-74. doi:10.1016/j.avsg.2023.05.046
6. Goueffic Y, Picquet J, Schneider F, et al. A Randomized Trial Comparing Polymer Versus Suture- Based Vascular Closure Devices for Arterial Closure Following Lower-Limb Arterial Endovascular Revascularization. *Cardiovasc Intervent Radiol*. Dec 2021;44(12):1883-1892. doi:10.1007/s00270-021-02940-z
7. Ha TH, Lee SH, Park SJ, Chung HH. Could real-time ultrasonography guidance be useful for the effective deployment of FemoSeal in common femoral arteriotomy? *Ultrasonography*. Jul 2021;40(3):449-454. doi:10.14366/usg.20127
8. Hassu J, Meyborg M, Gebauer K, Makowski L, Malyar N. 51. Embolization of closure device particle as a rare case of new-onset claudication - Normal diagnostic findings despite limiting symptoms. A Diagnostic Challenge. *Conference Abstract. Vasa*. Nov 2022;51(Suppl 109):1-21. doi:10.1024/0301-1526/a001037
9. Kaszczewski P, Kozubek H, Ostrowski T, et al. Rapid Formation and Hybrid Treatment of a Large Superior Mesenteric Artery Aneurysm. *Article. Am J Case Rep*. Oct 29 2023;24:e939558. doi:10.12659/AJCR.939558

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Funktionelle Tierversuche vor dem Inverkehrbringen	Produktversion
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vasculature System. Functional Implantation Test in Sheep 12-, 15-, and 18-month Time Periods. (R1051-1), 2000, Radi Medical System AB (performed by Biomatech, France).	11200 (Originalprodukt)
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep 12-, 24-, and 36-Week Time-Periods. (R1050-1), 2002, Radi Medical System AB (performed by Biomatech, France).	11200 (Originalprodukt)
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep and in the Pig for 5 Weeks, (R1052-01) 2002, Radi Medical System AB (performed by Biomatech, France).	11200 (Originalprodukt)
Erste klinische Humanstudie	
Multi-Centre Clinical Trial of the FemoSeal Vascular Closure System for Sealing Femoral Arterial Punctures After Diagnostic/Interventional Cardiology Procedures (EU-SE-001); 2004, Radi Medical Systems AB.	11200 (Originalprodukt)
Abgeschlossene klinische Studien nach dem Inverkehrbringen	
PMS FemoSeal® VCS #11202 (R1756-01) (2006)	C11202 (aktuelles Produkt)
Confirmatory Assessment of FemoSeal™ Vascular Closure System in 7F-Sheath. (2008)	C11202 (aktuelles Produkt)
Post-Marketing Surveillance (PMS) Study FemoSeal® VCS #11202 (R1756-03) (2009)	C11202 (aktuelles Produkt)
The FemoSeal Vascular Closure System (VCS) Registry: A Prospective, Multi-Center, Observational Study in Europe (2022)	C11202 (aktuelles Produkt)
Geplante/laufende klinische Studien nach dem Inverkehrbringen	
FemoSeal™ Vascular Closure System PMCF Clinical Survey Protocol: Addressing Gaps Identified in CER-FS-2023	C11202 (aktuelles Produkt)
Klinische Nachbeobachtungen nach dem Inverkehrbringen	
20210292 - FemoSeal Vascular Closure System Clinical Survey Report (2020)	C11202 (aktuelles Produkt)
20210298 – Guidewire Clinical Survey Report (2020)	C11202 (aktuelles Produkt)
FemoSeal™ Vascular Closure System PMCF Clinical Survey Report: Addressing Gaps Identified in CER-FS-2021 and CER-FS-2022 (2024)	C11202 (aktuelles Produkt)

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Pruebas funcionales con animales previas a la comercialización	Versión del dispositivo
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vasculature System. Functional Implantation Test in Sheep 12-, 15-, and 18-month Time Periods. (R1051-1), 2000, Radi Medical System AB (performed by Biomatech, France).	11200 (dispositivo original)
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep 12-, 24-, and 36-Week Time-Periods. (R1050-1), 2002, Radi Medical System AB (performed by Biomatech, France).	11200 (dispositivo original)
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep and in the Pig for 5 Weeks, (R1052-01) 2002, Radi Medical System AB (performed by Biomatech, France).	11200 (dispositivo original)
Primer estudio clínico en seres humanos	
Multi-Centre Clinical Trial of the FemoSeal Vascular Closure System for Sealing Femoral Arterial Punctures After Diagnostic/Interventional Cardiology Procedures (EU-SE-001); 2004, Radi Medical Systems AB.	11200 (dispositivo original)
Estudios clínicos posteriores a la comercialización finalizados	
PMS FemoSeal® VCS #11202 (R1756-01) (2006)	C11202 (dispositivo actual)
Confirmatory Assessment of FemoSeal™ Vascular Closure System in 7F-Sheath. (2008)	C11202 (dispositivo actual)
Post-Marketing Surveillance (PMS) Study FemoSeal® VCS #11202 (R1756-03) (2009)	C11202 (dispositivo actual)
The FemoSeal Vascular Closure System (VCS) Registry: A Prospective, Multi-Center, Observational Study in Europe (2022)	C11202 (dispositivo actual)
Estudios clínicos posteriores a la comercialización planificados/en curso	
FemoSeal™ Vascular Closure System PMCF Clinical Survey Protocol: Addressing Gaps Identified in CER-FS-2023	C11202 (dispositivo actual)

Actividades de seguimiento clínico posterior a la comercialización	
20210292 - FemoSeal Vascular Closure System Clinical Survey Report (2020)	C11202 (dispositivo actual)
20210298 – Guidewire Clinical Survey Report (2020)	C11202 (dispositivo actual)
FemoSeal™ Vascular Closure System PMCF Clinical Survey Report: Addressing Gaps Identified in CER-FS-2021 and CER-FS-2022 (2024)	C11202 (dispositivo actual)

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Turustamiseelised funktsionaalsed loomkatsed	Seadme versioon
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vasculature System. Functional Implantation Test in Sheep 12-, 15-, and 18-month Time Periods. (R1051-1), 2000, Radi Medical System AB (performed by Biomatech, France).	11200 (algne seade)
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep 12-, 24-, and 36-Week Time-Periods. (R1050-1), 2002, Radi Medical System AB (performed by Biomatech, France).	11200 (algne seade)
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep and in the Pig for 5 Weeks, (R1052-01) 2002, Radi Medical System AB (performed by Biomatech, France).	11200 (algne seade)
Esimene kliiniline uuring inimesel	
Multi-Centre Clinical Trial of the FemoSeal Vascular Closure System for Sealing Femoral Arterial Punctures After Diagnostic/Interventional Cardiology Procedures (EU-SE-001); 2004, Radi Medical Systems AB.	11200 (algne seade)
Lõpetatud turustamisjärgsed kliinilised uuringud	
PMS FemoSeal® VCS #11202 (R1756-01) (2006)	C11202 (praegune seade)
Confirmatory Assessment of FemoSeal™ Vascular Closure System in 7F-Sheath. (2008)	C11202 (praegune seade)
Post-Marketing Surveillance (PMS) Study FemoSeal® VCS #11202 (R1756-03) (2009)	C11202 (praegune seade)
The FemoSeal Vascular Closure System (VCS) Registry: A Prospective, Multi-Center, Observational Study in Europe (2022)	C11202 (praegune seade)
Plaanitavad/käimasolevad turustamisjärgsed kliinilised uuringud	
FemoSeal™ Vascular Closure System PMCF Clinical Survey Protocol: Addressing Gaps Identified in CER-FS-2023	C11202 (praegune seade)
Turustamisjärgse kliinilise järelkontrolli toimingud	
20210292 - FemoSeal Vascular Closure System Clinical Survey Report (2020)	C11202 (praegune seade)
20210298 – Guidewire Clinical Survey Report (2020)	C11202 (praegune seade)
FemoSeal™ Vascular Closure System PMCF Clinical Survey Report: Addressing Gaps Identified in CER-FS-2021 and CER-FS-2022 (2024)	C11202 (praegune seade)

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Toimintaa koskevat eläinkokeet ennen markkinoille saattamista	Laiteversio
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vasculature System. Functional Implantation Test in Sheep 12-, 15-, and 18-month Time Periods. (R1051-1), 2000, Radi Medical System AB (performed by Biomatech, France).	11200 (alkuperäinen laite)
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep 12-, 24-, and 36-Week Time-Periods. (R1050-1), 2002, Radi Medical System AB (performed by Biomatech, France).	11200 (alkuperäinen laite)
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep and in the Pig for 5 Weeks, (R1052-01) 2002, Radi Medical System AB (performed by Biomatech, France).	11200 (alkuperäinen laite)
Ensimmäinen kliininen ihmistutkimus	
Multi-Centre Clinical Trial of the FemoSeal Vascular Closure System for Sealing Femoral Arterial Punctures After Diagnostic/Interventional Cardiology Procedures (EU-SE-001); 2004, Radi Medical Systems AB.	11200 (alkuperäinen laite)
Markkinoille saattamisen jälkeen suoritettut kliiniset tutkimukset	
PMS FemoSeal® VCS #11202 (R1756-01) (2006)	C11202 (nykyinen laite)
Confirmatory Assessment of FemoSeal™ Vascular Closure System in 7F-Sheath. (2008)	C11202 (nykyinen laite)
Post-Marketing Surveillance (PMS) Study FemoSeal® VCS #11202 (R1756-03) (2009)	C11202 (nykyinen laite)
The FemoSeal Vascular Closure System (VCS) Registry: A Prospective, Multi-Center, Observational Study in Europe (2022)	C11202 (nykyinen laite)
Suunnitellut / meneillään olevat markkinoille saattamisen jälkeiset kliiniset tutkimukset	
FemoSeal™ Vascular Closure System PMCF Clinical Survey Protocol: Addressing Gaps Identified in CER-FS-2023	C11202 (nykyinen laite)

Kliiniset seurantatoimet markkinoille saattamisen jälkeen	
20210292 - FemoSeal Vascular Closure System Clinical Survey Report (2020)	C11202 (nykyinen laite)
20210298 – Guidewire Clinical Survey Report (2020)	C11202 (nykyinen laite)
FemoSeal™ Vascular Closure System PMCF Clinical Survey Report: Addressing Gaps Identified in CER-FS-2021 and CER-FS-2022 (2024)	C11202 (nykyinen laite)

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10. Kaszczewski P, Leszczynski J, Elwertowski M, Maciag R, Chudzinski W, Galazka Z. Combined Treatment of Multiple Splanchnic Artery Aneurysms Secondary to Median Arcuate Ligament Syndrome: A Case Study and Review of the Literature. Article. Am J Case Rep. Aug 11 2020;21:e926074. doi:10.12659/AJCR.926074
11. Kennedy SA, Rajan DK, Bassett P, Tan KT, Jaber A, Mafeld S. Complication rates associated with antegrade use of vascular closure devices: a systematic review and pooled analysis. J Vasc Surg. Feb 2021;73(2):722-730 e1. doi:10.1016/j.jvs.2020.08.133
12. Langouet Q, Martinez R, Saint-Etienne C, et al. Incidence, predictors, impact, and treatment of vascular complications after transcatheter aortic valve implantation in a modern prospective cohort under real conditions. Article. J Vasc Surg. Dec 2020;72(6):2120-2129 e2. doi:10.1016/j.jvs.2020.03.035
13. Mach M, Maciejewski K, Ostrowski T, Maciag R, Sajdek M, Galazka Z. Endovascular Treatment of a Bilateral, Ruptured Angiomyolipoma in a Patient With Tuberous Sclerosis Complex. Cureus. Aug 2024;16(8):e66200. doi:10.7759/cureus.66200
14. Mayer K, Gewalt S, Morath T, et al. Comparison of Vascular Closure Devices vs Manual Compression After Femoral Artery Puncture in Patients on Oral Anticoagulation - Post Hoc Analysis of the ISAR-CLOSURE Trial. The Journal of invasive cardiology. Sep 2021 2021;33(9):E709-E715.
15. Messana G, Ambrosi L, Moramarco LP, Cionfoli N, Maestri M, Quaretti P. Testicular artery originating from the inferior mesenteric artery: an alert for interventionalists – A case report. Article. Radiology Case Reports. 2021;16(9):2710-2713. doi:10.1016/j.radcr.2021.06.059
16. Moramarco LP, Capodaglio CA, Quaretti P, et al. Multivessel endovascular therapy for undiagnosed vascular type Ehlers-Danlos syndrome. Successful percutaneous transcatheter coil embolization of hepatic artery pseudoaneurysm with stenting of right renal and iliac arteries in emergency setting. BJR Case Rep. Dec 1 2020;6(4):20200025. doi:10.1259/bjrcr.20200025
17. Morosetti D, Chiocchi M, Argiro R, et al. Endovascular treatment of calcific lesions of the common femoral artery using atherectomy device associated with scoring balloon angioplasty in diabetic patients with high "major amputation" risk. Vascular. Jun 2022;30(3):463-473. doi:10.1177/17085381211019244
18. Noory E, Bohme T, Krause L, et al. Evaluation of the MYNX CONTROL Arterial Closure System for Achieving Primary Hemostasis after Arterial Femoral Access Following Peripheral Arterial Interventions, Compared to the FemoSeal(TM) Closure System. J Clin Med. Aug 12 2023;12(16):5255. doi:10.3390/jcm12165255

19. Oliny A, Bosse C, Charbonneau P, Haulon S, Fabre D. Closure Device Migration: An Unusual Cause of Acute Limb Ischaemia Following a Simple Endovascular Procedure. Article. EJVES Vasc Forum. 2024;61:113-115. doi:10.1016/j.ejvsf.2024.05.006
20. Pane F, Coppola M, Borzelli A, et al. A singular case of massive urethrorrhagia solved by transarterial embolization. Article. Egyptian Journal of Radiology and Nuclear Medicine. 2020;51(1):168. doi:10.1186/s43055-020-00299-8
21. Pang N, Gao J, Zhang B, et al. Vascular Closure Devices versus Manual Compression in Cardiac Interventional Procedures: Systematic Review and Meta-Analysis. Cardiovasc Ther. 2022 2022;2022:8569188. doi:10.1155/2022/8569188
22. Tagliaferro FB, Orgera G, Mascagni L, et al. FemoSeal((R)) vascular closure device for antegrade common femoral artery access: Safety and technical notes. J Vasc Access. Jan 2020;21(1):79- 85. doi:10.1177/1129729819854593
23. Tan R, Sieunarine K. Superficial Femoral Artery Pseudoaneurysm as a Delayed Complication of TurboHawk Atherectomy. Annals of Vascular Surgery. Oct 2020 2020;68:568.e17-568.e21. doi:10.1016/j.avsg.2020.04.006

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Essais fonctionnels pré-commercialisation sur des animaux	Version du dispositif
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vasculature System. Functional Implantation Test in Sheep 12-, 15-, and 18-month Time Periods. (R1051-1), 2000, Radi Medical System AB (performed by Biomatech, France).	11200 (dispositif d'origine)
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep 12-, 24-, and 36-Week Time-Periods. (R1050-1), 2002, Radi Medical System AB (performed by Biomatech, France).	11200 (dispositif d'origine)
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep and in the Pig for 5 Weeks, (R1052-01) 2002, Radi Medical System AB (performed by Biomatech, France).	11200 (dispositif d'origine)
Première étude clinique chez l'homme	
Multi-Centre Clinical Trial of the FemoSeal Vascular Closure System for Sealing Femoral Arterial Punctures After Diagnostic/Interventional Cardiology Procedures (EU-SE-001); 2004, Radi Medical Systems AB.	11200 (dispositif d'origine)
Études cliniques post-commercialisation terminées	
PMS FemoSeal® VCS #11202 (R1756-01) (2006)	C11202 (appareil actuel)
Confirmatory Assessment of FemoSeal™ Vascular Closure System in 7F-Sheath. (2008)	C11202 (appareil actuel)
Post-Marketing Surveillance (PMS) Study FemoSeal® VCS #11202 (R1756-03) (2009)	C11202 (appareil actuel)
The FemoSeal Vascular Closure System (VCS) Registry: A Prospective, Multi-Center, Observational Study in Europe (2022)	C11202 (appareil actuel)
Études cliniques post-commercialisation planifiées/en cours	
FemoSeal™ Vascular Closure System PMCF Clinical Survey Protocol: Addressing Gaps Identified in CER-FS-2023	C11202 (appareil actuel)

Activités de suivi clinique post-commercialisation	
20210292 - FemoSeal Vascular Closure System Clinical Survey Report (2020)	C11202 (appareil actuel)
20210298 – Guidewire Clinical Survey Report (2020)	C11202 (appareil actuel)
FemoSeal™ Vascular Closure System PMCF Clinical Survey Report: Addressing Gaps Identified in CER-FS-2021 and CER-FS-2022 (2024)	C11202 (appareil actuel)

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2.0 BIBLIOGRAPHIE

1. Alabduljabbar M, Alhamdan B, Alabdrabalnabi M. Access site complications of manual compression versus closure devices after lower limb revascularization: a systematic review. *International Journal of Medicine in Developing Countries*. 2022;6(12):1572-1580. doi:10.24911/ijmdc.51-1666716459
2. Chugh Y, Bavishi C, Mojadidi MK, et al. Safety of transradial access compared to transfemoral access with hemostatic devices (vessel plugs and suture devices) after percutaneous coronary interventions: A systematic review and meta-analysis. *Catheter Cardiovasc Interv*. Aug 2020;96(2):285-295. doi:10.1002/ccd.29061
3. Ciprian Cacuci A, Krankenberg H, Ingwersen M, et al. Access Site Complications of Peripheral Endovascular Procedures: A Large, Prospective Registry on Predictors and Consequences. *J Endovasc Ther*. Oct 2021;28(5):746-754. doi:10.1177/15266028211025044
4. Faraj J, Daneshmand A, Mwiapatayi BP. Large spontaneous median sacral artery pseudoaneurysm masquerading as hip pain. *Article. J Vasc Surg Cases Innov Tech*. Jun 2020;6(2):224-227. doi:10.1016/j.jvscit.2020.02.012
5. Gabrielli R, Siani A, Smedile G, Rizzo AR, Accrocca F, Bartoli S. Carotid Artery Stenting versus Carotid Endarterectomy in Terms of Neuroprotection DW-MRI Detected and Neuropsychological Assessment Impairment. *Article. Ann Vasc Surg*. Jan 2024;98:68-74. doi:10.1016/j.avsg.2023.05.046
6. Goueffic Y, Picquet J, Schneider F, et al. A Randomized Trial Comparing Polymer Versus Suture- Based Vascular Closure Devices for Arterial Closure Following Lower-Limb Arterial Endovascular Revascularization. *Cardiovasc Intervent Radiol*. Dec 2021;44(12):1883-1892. doi:10.1007/s00270-021-02940-z
7. Ha TH, Lee SH, Park SJ, Chung HH. Could real-time ultrasonography guidance be useful for the effective deployment of FemoSeal in common femoral arteriotomy? *Ultrasonography*. Jul 2021;40(3):449-454. doi:10.14366/usg.20127
8. Hassu J, Meyborg M, Gebauer K, Makowski L, Malyar N. 51. Embolization of closure device particle as a rare case of new-onset claudication - Normal diagnostic findings despite limiting symptoms. A Diagnostic Challenge. *Conference Abstract. Vasa*. Nov 2022;51(Suppl 109):1-21. doi:10.1024/0301-1526/a001037
9. Kaszczewski P, Kozubek H, Ostrowski T, et al. Rapid Formation and Hybrid Treatment of a Large Superior Mesenteric Artery Aneurysm. *Article. Am J Case Rep*. Oct 29 2023;24:e939558. doi:10.12659/AJCR.939558

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1. Alabduljabbar M, Alhamdan B, Alabdrabalnabi M. Access site complications of manual compression versus closure devices after lower limb revascularization: a systematic review. *International Journal of Medicine in Developing Countries*. 2022;6(12):1572-1580. doi:10.24911/ijmdc.51-1666716459
2. Chugh Y, Bavishi C, Mojadidi MK, et al. Safety of transradial access compared to transfemoral access with hemostatic devices (vessel plugs and suture devices) after percutaneous coronary interventions: A systematic review and meta-analysis. *Catheter Cardiovasc Interv*. Aug 2020;96(2):285-295. doi:10.1002/ccd.29061
3. Ciprian Cacuci A, Krankenberg H, Ingwersen M, et al. Access Site Complications of Peripheral Endovascular Procedures: A Large, Prospective Registry on Predictors and Consequences. *J Endovasc Ther*. Oct 2021;28(5):746-754. doi:10.1177/15266028211025044
4. Faraj J, Daneshmand A, Mwipatayi BP. Large spontaneous median sacral artery pseudoaneurysm masquerading as hip pain. *Article. J Vasc Surg Cases Innov Tech*. Jun 2020;6(2):224-227. doi:10.1016/j.jvscit.2020.02.012
5. Gabrielli R, Siani A, Smedile G, Rizzo AR, Accrocca F, Bartoli S. Carotid Artery Stenting versus Carotid Endarterectomy in Terms of Neuroprotection DW-MRI Detected and Neuropsychological Assessment Impairment. *Article. Ann Vasc Surg*. Jan 2024;98:68-74. doi:10.1016/j.avsg.2023.05.046
6. Goueffic Y, Picquet J, Schneider F, et al. A Randomized Trial Comparing Polymer Versus Suture- Based Vascular Closure Devices for Arterial Closure Following Lower-Limb Arterial Endovascular Revascularization. *Cardiovasc Intervent Radiol*. Dec 2021;44(12):1883-1892. doi:10.1007/s00270-021-02940-z
7. Ha TH, Lee SH, Park SJ, Chung HH. Could real-time ultrasonography guidance be useful for the effective deployment of FemoSeal in common femoral arteriotomy? *Ultrasonography*. Jul 2021;40(3):449-454. doi:10.14366/usg.20127
8. Hassu J, Meyborg M, Gebauer K, Makowski L, Malyar N. 51. Embolization of closure device particle as a rare case of new-onset claudication - Normal diagnostic findings despite limiting symptoms. A Diagnostic Challenge. *Conference Abstract. Vasa*. Nov 2022;51(Suppl 109):1-21. doi:10.1024/0301-1526/a001037
9. Kaszczewski P, Kozubek H, Ostrowski T, et al. Rapid Formation and Hybrid Treatment of a Large Superior Mesenteric Artery Aneurysm. *Article. Am J Case Rep*. Oct 29 2023;24:e939558. doi:10.12659/AJCR.939558

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Tablica 2.3 Sažetak ključnih ispitivanja

Funkcionalna ispitivanja na životinjama prije stavljanja u promet	Verzija proizvoda
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vasculature System. Functional Implantation Test in Sheep 12-, 15-, and 18-month Time Periods. (R1051-1), 2000, Radi Medical System AB (performed by Biomatech, France).	11200 (originalni proizvod)
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep 12-, 24-, and 36-Week Time-Periods. (R1050-1), 2002, Radi Medical System AB (performed by Biomatech, France).	11200 (originalni proizvod)
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep and in the Pig for 5 Weeks, (R1052-01) 2002, Radi Medical System AB (performed by Biomatech, France).	11200 (originalni proizvod)
Prvo kliničko ispitivanje na ljudima	
Multi-Centre Clinical Trial of the FemoSeal Vascular Closure System for Sealing Femoral Arterial Punctures After Diagnostic/Interventional Cardiology Procedures (EU-SE-001); 2004, Radi Medical Systems AB.	11200 (originalni proizvod)
Završena klinička ispitivanja nakon stavljanja u promet	
PMS FemoSeal® VCS #11202 (R1756-01) (2006)	C11202 (trenutačni proizvod)
Confirmatory Assessment of FemoSeal™ Vascular Closure System in 7F-Sheath. (2008)	C11202 (trenutačni proizvod)
Post-Marketing Surveillance (PMS) Study FemoSeal® VCS #11202 (R1756-03) (2009)	C11202 (trenutačni proizvod)
The FemoSeal Vascular Closure System (VCS) Registry: A Prospective, Multi-Center, Observational Study in Europe (2022)	C11202 (trenutačni proizvod)
Planirana/tekuća klinička ispitivanja nakon stavljanja u promet	
FemoSeal™ Vascular Closure System PMCF Clinical Survey Protocol: Addressing Gaps Identified in CER-FS-2023	C11202 (trenutačni proizvod)
Aktivnosti kliničkog praćenja nakon stavljanja u promet	
20210292 - FemoSeal Vascular Closure System Clinical Survey Report (2020)	C11202 (trenutačni proizvod)
20210298 – Guidewire Clinical Survey Report (2020)	C11202 (trenutačni proizvod)
FemoSeal™ Vascular Closure System PMCF Clinical Survey Report: Addressing Gaps Identified in CER-FS-2021 and CER-FS-2022 (2024)	C11202 (trenutačni proizvod)

Tablica 2.4 u nastavku sadrži opis svih ispitivanja sustava FemoSeal™ VCS prije stavljanja u promet.

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10. Kaszczewski P, Leszczynski J, Elwertowski M, Maciag R, Chudzinski W, Galazka Z. Combined Treatment of Multiple Splanchnic Artery Aneurysms Secondary to Median Arcuate Ligament Syndrome: A Case Study and Review of the Literature. Article. Am J Case Rep. Aug 11 2020;21:e926074. doi:10.12659/AJCR.926074
11. Kennedy SA, Rajan DK, Bassett P, Tan KT, Jaber A, Mafeld S. Complication rates associated with antegrade use of vascular closure devices: a systematic review and pooled analysis. J Vasc Surg. Feb 2021;73(2):722-730 e1. doi:10.1016/j.jvs.2020.08.133
12. Langouet Q, Martinez R, Saint-Etienne C, et al. Incidence, predictors, impact, and treatment of vascular complications after transcatheter aortic valve implantation in a modern prospective cohort under real conditions. Article. J Vasc Surg. Dec 2020;72(6):2120-2129 e2. doi:10.1016/j.jvs.2020.03.035
13. Mach M, Maciejewski K, Ostrowski T, Maciag R, Sajdek M, Galazka Z. Endovascular Treatment of a Bilateral, Ruptured Angiomyolipoma in a Patient With Tuberos Sclerosis Complex. Cureus. Aug 2024;16(8):e66200. doi:10.7759/cureus.66200
14. Mayer K, Gewalt S, Morath T, et al. Comparison of Vascular Closure Devices vs Manual Compression After Femoral Artery Puncture in Patients on Oral Anticoagulation - Post Hoc Analysis of the ISAR-CLOSURE Trial. The Journal of invasive cardiology. Sep 2021 2021;33(9):E709-E715.
15. Messana G, Ambrosi L, Moramarco LP, Cionfoli N, Maestri M, Quaretti P. Testicular artery originating from the inferior mesenteric artery: an alert for interventionalists – A case report. Article. Radiology Case Reports. 2021;16(9):2710-2713. doi:10.1016/j.radcr.2021.06.059
16. Moramarco LP, Capodaglio CA, Quaretti P, et al. Multivessel endovascular therapy for undiagnosed vascular type Ehlers-Danlos syndrome. Successful percutaneous transcatheter coil embolization of hepatic artery pseudoaneurysm with stenting of right renal and iliac arteries in emergency setting. BJR Case Rep. Dec 1 2020;6(4):20200025. doi:10.1259/bjrcr.20200025
17. Morosetti D, Chiocchi M, Argiro R, et al. Endovascular treatment of calcific lesions of the common femoral artery using atherectomy device associated with scoring balloon angioplasty in diabetic patients with high "major amputation" risk. Vascular. Jun 2022;30(3):463-473. doi:10.1177/17085381211019244
18. Noory E, Bohme T, Krause L, et al. Evaluation of the MYNX CONTROL Arterial Closure System for Achieving Primary Hemostasis after Arterial Femoral Access Following Peripheral Arterial Interventions, Compared to the FemoSeal(TM) Closure System. J Clin Med. Aug 12 2023;12(16):5255. doi:10.3390/jcm12165255

19. Oliny A, Bosse C, Charbonneau P, Haulon S, Fabre D. Closure Device Migration: An Unusual Cause of Acute Limb Ischaemia Following a Simple Endovascular Procedure. Article. EJVES Vasc Forum. 2024;61:113-115. doi:10.1016/j.ejvsf.2024.05.006
20. Pane F, Coppola M, Borzelli A, et al. A singular case of massive urethrorrhagia solved by transarterial embolization. Article. Egyptian Journal of Radiology and Nuclear Medicine. 2020;51(1):168. doi:10.1186/s43055-020-00299-8
21. Pang N, Gao J, Zhang B, et al. Vascular Closure Devices versus Manual Compression in Cardiac Interventional Procedures: Systematic Review and Meta-Analysis. Cardiovasc Ther. 2022 2022;2022:8569188. doi:10.1155/2022/8569188
22. Tagliaferro FB, Orgera G, Mascagni L, et al. FemoSeal((R)) vascular closure device for antegrade common femoral artery access: Safety and technical notes. J Vasc Access. Jan 2020;21(1):79- 85. doi:10.1177/1129729819854593
23. Tan R, Sieunarine K. Superficial Femoral Artery Pseudoaneurysm as a Delayed Complication of TurboHawk Atherectomy. Annals of Vascular Surgery. Oct 2020 2020;68:568.e17-568.e21. doi:10.1016/j.avsg.2020.04.006

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Test funzionali pre-commercializzazione sugli animali	Versione dispositivo
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vasculature System. Functional Implantation Test in Sheep 12-, 15-, and 18-month Time Periods. (R1051-1), 2000, Radi Medical System AB (performed by Biomatech, France).	11200 (dispositivo originale)
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep 12-, 24-, and 36-Week Time-Periods. (R1050-1), 2002, Radi Medical System AB (performed by Biomatech, France).	11200 (dispositivo originale)
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep and in the Pig for 5 Weeks, (R1052-01) 2002, Radi Medical System AB (performed by Biomatech, France).	11200 (dispositivo originale)
Primo studio clinico sull'uomo	
Multi-Centre Clinical Trial of the FemoSeal Vascular Closure System for Sealing Femoral Arterial Punctures After Diagnostic/Interventional Cardiology Procedures (EU-SE-001); 2004, Radi Medical Systems AB.	11200 (dispositivo originale)
Studi clinici post-commercializzazione completati	
PMS FemoSeal® VCS #11202 (R1756-01) (2006)	C11202 (dispositivo attuale)
Confirmatory Assessment of FemoSeal™ Vascular Closure System in 7F-Sheath. (2008)	C11202 (dispositivo attuale)
Post-Marketing Surveillance (PMS) Study FemoSeal® VCS #11202 (R1756-03) (2009)	C11202 (dispositivo attuale)
The FemoSeal Vascular Closure System (VCS) Registry: A Prospective, Multi-Center, Observational Study in Europe (2022)	C11202 (dispositivo attuale)
Studi clinici post-commercializzazione pianificati/in corso	
FemoSeal™ Vascular Closure System PMCF Clinical Survey Protocol: Addressing Gaps Identified in CER-FS-2023	C11202 (dispositivo attuale)

Attività di follow-up clinico post-commercializzazione	
20210292 - FemoSeal Vascular Closure System Clinical Survey Report (2020)	C11202 (dispositivo attuale)
20210298 – Guidewire Clinical Survey Report (2020)	C11202 (dispositivo attuale)
FemoSeal™ Vascular Closure System PMCF Clinical Survey Report: Addressing Gaps Identified in CER-FS-2021 and CER-FS-2022 (2024)	C11202 (dispositivo attuale)

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시판 전 기능성 동물 실험	장치 버전
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vasculature System. Functional Implantation Test in Sheep 12-, 15-, and 18-month Time Periods. (R1051-1), 2000, Radi Medical System AB (performed by Biomatech, France).	11200 (원래 장치)
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep 12-, 24-, and 36-Week Time-Periods. (R1050-1), 2002, Radi Medical System AB (performed by Biomatech, France).	11200 (원래 장치)
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep and in the Pig for 5 Weeks, (R1052-01) 2002, Radi Medical System AB (performed by Biomatech, France).	11200 (원래 장치)
최초의 사람 대상 임상 시험	
Multi-Centre Clinical Trial of the FemoSeal Vascular Closure System for Sealing Femoral Arterial Punctures After Diagnostic/Interventional Cardiology Procedures (EU-SE-001); 2004, Radi Medical Systems AB.	11200 (원래 장치)
시판 후 임상 연구 완료	
PMS FemoSeal® VCS #11202 (R1756-01) (2006)	C11202 (현재 장치)
Confirmatory Assessment of FemoSeal™ Vascular Closure System in 7F-Sheath. (2008)	C11202 (현재 장치)
Post-Marketing Surveillance (PMS) Study FemoSeal® VCS #11202 (R1756-03) (2009)	C11202 (현재 장치)
The FemoSeal Vascular Closure System (VCS) Registry: A Prospective, Multi-Center, Observational Study in Europe (2022)	C11202 (현재 장치)
계획/진행 중인 시판 후 임상시험	
FemoSeal™ Vascular Closure System PMCF Clinical Survey Protocol: Addressing Gaps Identified in CER-FS-2023	C11202 (현재 장치)
시판 후 임상 추적관찰 활동	
20210292 - FemoSeal Vascular Closure System Clinical Survey Report (2020)	C11202 (현재 장치)
20210298 – Guidewire Clinical Survey Report (2020)	C11202 (현재 장치)
FemoSeal™ Vascular Closure System PMCF Clinical Survey Report: Addressing Gaps Identified in CER-FS-2021 and CER-FS-2022 (2024)	C11202 (현재 장치)

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10. Kaszczewski P, Leszczynski J, Elwertowski M, Maciag R, Chudzinski W, Galazka Z. Combined Treatment of Multiple Splanchnic Artery Aneurysms Secondary to Median Arcuate Ligament Syndrome: A Case Study and Review of the Literature. Article. Am J Case Rep. Aug 11 2020;21:e926074. doi:10.12659/AJCR.926074
11. Kennedy SA, Rajan DK, Bassett P, Tan KT, Jaber A, Mafeld S. Complication rates associated with antegrade use of vascular closure devices: a systematic review and pooled analysis. J Vasc Surg. Feb 2021;73(2):722-730 e1. doi:10.1016/j.jvs.2020.08.133
12. Langouet Q, Martinez R, Saint-Etienne C, et al. Incidence, predictors, impact, and treatment of vascular complications after transcatheter aortic valve implantation in a modern prospective cohort under real conditions. Article. J Vasc Surg. Dec 2020;72(6):2120-2129 e2. doi:10.1016/j.jvs.2020.03.035
13. Mach M, Maciejewski K, Ostrowski T, Maciag R, Sajdek M, Galazka Z. Endovascular Treatment of a Bilateral, Ruptured Angiomyolipoma in a Patient With Tuberous Sclerosis Complex. Cureus. Aug 2024;16(8):e66200. doi:10.7759/cureus.66200
14. Mayer K, Gewalt S, Morath T, et al. Comparison of Vascular Closure Devices vs Manual Compression After Femoral Artery Puncture in Patients on Oral Anticoagulation - Post Hoc Analysis of the ISAR-CLOSURE Trial. The Journal of invasive cardiology. Sep 2021 2021;33(9):E709-E715.
15. Messana G, Ambrosi L, Moramarco LP, Cionfoli N, Maestri M, Quaretti P. Testicular artery originating from the inferior mesenteric artery: an alert for interventionalists – A case report. Article. Radiology Case Reports. 2021;16(9):2710-2713. doi:10.1016/j.radcr.2021.06.059
16. Moramarco LP, Capodaglio CA, Quaretti P, et al. Multivessel endovascular therapy for undiagnosed vascular type Ehlers-Danlos syndrome. Successful percutaneous transcatheter coil embolization of hepatic artery pseudoaneurysm with stenting of right renal and iliac arteries in emergency setting. BJR Case Rep. Dec 1 2020;6(4):20200025. doi:10.1259/bjrcr.20200025
17. Morosetti D, Chicchi M, Argiro R, et al. Endovascular treatment of calcific lesions of the common femoral artery using atherectomy device associated with scoring balloon angioplasty in diabetic patients with high "major amputation" risk. Vascular. Jun 2022;30(3):463-473. doi:10.1177/17085381211019244
18. Noory E, Bohme T, Krause L, et al. Evaluation of the MYNX CONTROL Arterial Closure System for Achieving Primary Hemostasis after Arterial Femoral Access Following Peripheral Arterial Interventions, Compared to the FemoSeal(TM) Closure System. J Clin Med. Aug 12 2023;12(16):5255. doi:10.3390/jcm12165255

19. Oliny A, Bosse C, Charbonneau P, Haulon S, Fabre D. Closure Device Migration: An Unusual Cause of Acute Limb Ischaemia Following a Simple Endovascular Procedure. Article. EJVES Vasc Forum. 2024;61:113-115. doi:10.1016/j.ejvsf.2024.05.006
20. Pane F, Coppola M, Borzelli A, et al. A singular case of massive urethrorrhagia solved by transarterial embolization. Article. Egyptian Journal of Radiology and Nuclear Medicine. 2020;51(1):168. doi:10.1186/s43055-020-00299-8
21. Pang N, Gao J, Zhang B, et al. Vascular Closure Devices versus Manual Compression in Cardiac Interventional Procedures: Systematic Review and Meta-Analysis. Cardiovasc Ther. 2022 2022;2022:8569188. doi:10.1155/2022/8569188
22. Tagliaferro FB, Orgera G, Mascagni L, et al. FemoSeal((R)) vascular closure device for antegrade common femoral artery access: Safety and technical notes. J Vasc Access. Jan 2020;21(1):79- 85. doi:10.1177/1129729819854593
23. Tan R, Sieunarine K. Superficial Femoral Artery Pseudoaneurysm as a Delayed Complication of TurboHawk Atherectomy. Annals of Vascular Surgery. Oct 2020 2020;68:568.e17-568.e21. doi:10.1016/j.avsg.2020.04.006

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Pirmspārdošanas funkcionālā testēšana ar dzīvniekiem	Ierīces versija
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vasculature System. Functional Implantation Test in Sheep 12-, 15-, and 18-month Time Periods. (R1051-1), 2000, Radi Medical System AB (performed by Biomatech, France).	11200 (oriģinālā ierīce)
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep 12-, 24-, and 36-Week Time-Periods. (R1050-1), 2002, Radi Medical System AB (performed by Biomatech, France).	11200 (oriģinālā ierīce)
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep and in the Pig for 5 Weeks, (R1052-01) 2002, Radi Medical System AB (performed by Biomatech, France).	11200 (oriģinālā ierīce)
Pirmais klīniskais pētījums ar cilvēkiem	
Multi-Centre Clinical Trial of the FemoSeal Vascular Closure System for Sealing Femoral Arterial Punctures After Diagnostic/Interventional Cardiology Procedures (EU-SE-001); 2004, Radi Medical Systems AB.	11200 (oriģinālā ierīce)
Pabeigti pēcpārdošanas klīniskie pētījumi	
PMS FemoSeal® VCS #11202 (R1756-01) (2006)	C11202 (pašreizējā ierīce)
Confirmatory Assessment of FemoSeal™ Vascular Closure System in 7F-Sheath. (2008)	C11202 (pašreizējā ierīce)
Post-Marketing Surveillance (PMS) Study FemoSeal® VCS #11202 (R1756-03) (2009)	C11202 (pašreizējā ierīce)
The FemoSeal Vascular Closure System (VCS) Registry: A Prospective, Multi-Center, Observational Study in Europe (2022)	C11202 (pašreizējā ierīce)
Plānotie/notiekošie pēcpārdošanas klīniskie pētījumi	
FemoSeal™ Vascular Closure System PMCF Clinical Survey Protocol: Addressing Gaps Identified in CER-FS-2023	C11202 (pašreizējā ierīce)
Pēcpārdošanas klīniskās novērošanas pasākumi	
20210292 - FemoSeal Vascular Closure System Clinical Survey Report (2020)	C11202 (pašreizējā ierīce)
20210298 – Guidewire Clinical Survey Report (2020)	C11202 (pašreizējā ierīce)
FemoSeal™ Vascular Closure System PMCF Clinical Survey Report: Addressing Gaps Identified in CER-FS-2021 and CER-FS-2022 (2024)	C11202 (pašreizējā ierīce)

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Funkcinis bandymas su gyvūnais prieš pateikiant rinkai	Įtaiso versija
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vasculature System. Functional Implantation Test in Sheep 12-, 15-, and 18-month Time Periods. (R1051-1), 2000, Radi Medical System AB (performed by Biomatech, France).	11200 (originalus įtaisas)
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep 12-, 24-, and 36-Week Time-Periods. (R1050-1), 2002, Radi Medical System AB (performed by Biomatech, France).	11200 (originalus įtaisas)
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep and in the Pig for 5 Weeks, (R1052-01) 2002, Radi Medical System AB (performed by Biomatech, France).	11200 (originalus įtaisas)
Pirmasis klinikinis tyrimas su žmonėmis	
Multi-Centre Clinical Trial of the FemoSeal Vascular Closure System for Sealing Femoral Arterial Punctures After Diagnostic/Interventional Cardiology Procedures (EU-SE-001); 2004, Radi Medical Systems AB.	11200 (originalus įtaisas)
Po pateikimo į rinką atlikti klinikiniai tyrimai	
PMS FemoSeal® VCS #11202 (R1756-01) (2006)	C11202 (dabartinis įtaisas)
Confirmatory Assessment of FemoSeal™ Vascular Closure System in 7F-Sheath. (2008)	C11202 (esamas prietaisas)
Post-Marketing Surveillance (PMS) Study FemoSeal® VCS #11202 (R1756-03) (2009)	C11202 (esamas prietaisas)
The FemoSeal Vascular Closure System (VCS) Registry: A Prospective, Multi-Center, Observational Study in Europe (2022)	C11202 (esamas prietaisas)
Planuojami / vykdomi klinikiniai tyrimai po pateikimo rinkai	
FemoSeal™ Vascular Closure System PMCF Clinical Survey Protocol: Addressing Gaps Identified in CER-FS-2023	C11202 (esamas prietaisas)
Tolesnė klinikinė veikla po pateikimo rinkai	
20210292 - FemoSeal Vascular Closure System Clinical Survey Report (2020)	C11202 (esamas prietaisas)

20210298 – Guidewire Clinical Survey Report (2020)	C11202 (esamas prietaisas)
FemoSeal™ Vascular Closure System PMCF Clinical Survey Report: Addressing Gaps Identified in CER-FS-2021 and CER-FS-2022 (2024)	C11202 (esamas prietaisas)

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2.0 BIBLIOGRAFIJA

1. Alabduljabbar M, Alhamdan B, Alabdrabalnabi M. Access site complications of manual compression versus closure devices after lower limb revascularization: a systematic review. *International Journal of Medicine in Developing Countries*. 2022;6(12):1572-1580. doi:10.24911/ijmdc.51-1666716459
2. Chugh Y, Bavishi C, Mojadidi MK, et al. Safety of transradial access compared to transfemoral access with hemostatic devices (vessel plugs and suture devices) after percutaneous coronary interventions: A systematic review and meta-analysis. *Catheter Cardiovasc Interv*. Aug 2020;96(2):285-295. doi:10.1002/ccd.29061
3. Ciprian Cacuci A, Krankenberg H, Ingwersen M, et al. Access Site Complications of Peripheral Endovascular Procedures: A Large, Prospective Registry on Predictors and Consequences. *J Endovasc Ther*. Oct 2021;28(5):746-754. doi:10.1177/15266028211025044
4. Faraj J, Daneshmand A, Mwiratayi BP. Large spontaneous median sacral artery pseudoaneurysm masquerading as hip pain. Article. *J Vasc Surg Cases Innov Tech*. Jun 2020;6(2):224-227. doi:10.1016/j.jvscit.2020.02.012
5. Gabrielli R, Siani A, Smedile G, Rizzo AR, Accrocca F, Bartoli S. Carotid Artery Stenting versus Carotid Endarterectomy in Terms of Neuroprotection DW-MRI Detected and Neuropsychological Assessment Impairment. Article. *Ann Vasc Surg*. Jan 2024;98:68-74. doi:10.1016/j.avsg.2023.05.046
6. Goueffic Y, Picquet J, Schneider F, et al. A Randomized Trial Comparing Polymer Versus Suture- Based Vascular Closure Devices for Arterial Closure Following Lower-Limb Arterial Endovascular Revascularization. *Cardiovasc Intervent Radiol*. Dec 2021;44(12):1883-1892. doi:10.1007/s00270-021-02940-z
7. Ha TH, Lee SH, Park SJ, Chung HH. Could real-time ultrasonography guidance be useful for the effective deployment of FemoSeal in common femoral arteriotomy? *Ultrasonography*. Jul 2021;40(3):449-454. doi:10.14366/usg.20127
8. Hassu J, Meyborg M, Gebauer K, Makowski L, Malyar N. 51. Embolization of closure device particle as a rare case of new-onset claudication - Normal diagnostic findings despite limiting symptoms. A Diagnostic Challenge. Conference Abstract. *Vasa*. Nov 2022;51(Suppl 109):1-21. doi:10.1024/0301-1526/a001037
9. Kaszczewski P, Kozubek H, Ostrowski T, et al. Rapid Formation and Hybrid Treatment of a Large Superior Mesenteric Artery Aneurysm. Article. *Am J Case Rep*. Oct 29 2023;24:e939558. doi:10.12659/AJCR.939558

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Test funkčnosti na zvieratách pred uvedením na trh	Verzia pomôcky
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vasculature System. Functional Implantation Test in Sheep 12-, 15-, and 18-month Time Periods. (R1051-1), 2000, Radi Medical System AB (performed by Biomatech, France).	11200 (pôvodná pomôcka)
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep 12-, 24-, and 36-Week Time-Periods. (R1050-1), 2002, Radi Medical System AB (performed by Biomatech, France).	11200 (pôvodná pomôcka)
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep and in the Pig for 5 Weeks, (R1052-01) 2002, Radi Medical System AB (performed by Biomatech, France).	11200 (pôvodná pomôcka)
Prvá klinická štúdia na ľuďoch	
Multi-Centre Clinical Trial of the FemoSeal Vascular Closure System for Sealing Femoral Arterial Punctures After Diagnostic/Interventional Cardiology Procedures (EU-SE-001); 2004, Radi Medical Systems AB.	11200 (pôvodná pomôcka)
Dokončené klinické štúdie po uvedení na trh	
PMS FemoSeal® VCS #11202 (R1756-01) (2006)	C11202 (aktuálna pomôcka)
Confirmatory Assessment of FemoSeal™ Vascular Closure System in 7F-Sheath. (2008)	C11202 (aktuálna pomôcka)
Post-Marketing Surveillance (PMS) Study FemoSeal® VCS #11202 (R1756-03) (2009)	C11202 (aktuálna pomôcka)
The FemoSeal Vascular Closure System (VCS) Registry: A Prospective, Multi-Center, Observational Study in Europe (2022)	C11202 (aktuálna pomôcka)
Plánované/prebiehajúce klinické štúdie po uvedení na trh	
FemoSeal™ Vascular Closure System PMCF Clinical Survey Protocol: Addressing Gaps Identified in CER-FS-2023	C11202 (aktuálna pomôcka)
Činnosti klinického sledovania po uvedení na trh	
20210292 - FemoSeal Vascular Closure System Clinical Survey Report (2020)	C11202 (aktuálna pomôcka)
20210298 – Guidewire Clinical Survey Report (2020)	C11202 (aktuálna pomôcka)
FemoSeal™ Vascular Closure System PMCF Clinical Survey Report: Addressing Gaps Identified in CER-FS-2021 and CER-FS-2022 (2024)	C11202 (aktuálna pomôcka)

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Funkcionalno testiranje na živalih pred dajanjem na trg	Različica pripomočka
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vasculature System. Functional Implantation Test in Sheep 12-, 15-, and 18-month Time Periods. (R1051-1), 2000, Radi Medical System AB (performed by Biomatech, France).	11200 (originalni pripomoček)
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep 12-, 24-, and 36-Week Time-Periods. (R1050-1), 2002, Radi Medical System AB (performed by Biomatech, France).	11200 (originalni pripomoček)
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep and in the Pig for 5 Weeks, (R1052-01) 2002, Radi Medical System AB (performed by Biomatech, France).	11200 (originalni pripomoček)
Prva klinična študija pri ljudeh	
Multi-Centre Clinical Trial of the FemoSeal Vascular Closure System for Sealing Femoral Arterial Punctures After Diagnostic/Interventional Cardiology Procedures (EU-SE-001); 2004, Radi Medical Systems AB.	11200 (originalni pripomoček)
Zaključene klinične študije po dajanju na trg	
PMS FemoSeal® VCS #11202 (R1756-01) (2006)	C11202 (trenutni pripomoček)
Confirmatory Assessment of FemoSeal™ Vascular Closure System in 7F-Sheath. (2008)	C11202 (trenutni pripomoček)
Post-Marketing Surveillance (PMS) Study FemoSeal® VCS #11202 (R1756-03) (2009)	C11202 (trenutni pripomoček)
The FemoSeal Vascular Closure System (VCS) Registry: A Prospective, Multi-Center, Observational Study in Europe (2022)	C11202 (trenutni pripomoček)
Načrtovane/potekajoče klinične študije po dajanju na trg	
FemoSeal™ Vascular Closure System PMCF Clinical Survey Protocol: Addressing Gaps Identified in CER-FS-2023	C11202 (trenutni pripomoček)
Dejavnosti kliničnega spremljanja po dajanju na trg	
20210292 - FemoSeal Vascular Closure System Clinical Survey Report (2020)	C11202 (trenutni pripomoček)

20210298 – Guidewire Clinical Survey Report (2020)	C11202 (trenutni pripomoček)
FemoSeal™ Vascular Closure System PMCF Clinical Survey Report: Addressing Gaps Identified in CER-FS-2021 and CER-FS-2022 (2024)	C11202 (trenutni pripomoček)

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Tabela 2.3 Rezime pivotalnih studija

Funkcionalno testiranje na životinjama pre plasiranja na tržište	Verzija sredstva
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vasculature System. Functional Implantation Test in Sheep 12-, 15-, and 18-month Time Periods. (R1051-1), 2000, Radi Medical System AB (performed by Biomatech, France).	11200 (originalno sredstvo)
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep 12-, 24-, and 36-Week Time-Periods. (R1050-1), 2002, Radi Medical System AB (performed by Biomatech, France).	11200 (originalno sredstvo)
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep and in the Pig for 5 Weeks, (R1052-01) 2002, Radi Medical System AB (performed by Biomatech, France).	11200 (originalno sredstvo)
Prva klinička studija na ljudima	
Multi-Centre Clinical Trial of the FemoSeal Vascular Closure System for Sealing Femoral Arterial Punctures After Diagnostic/Interventional Cardiology Procedures (EU-SE-001); 2004, Radi Medical Systems AB.	11200 (originalno sredstvo)
Završene kliničke studije nakon plasiranja na tržište	
PMS FemoSeal® VCS #11202 (R1756-01) (2006)	C11202 (trenutno sredstvo)
Confirmatory Assessment of FemoSeal™ Vascular Closure System in 7F-Sheath. (2008)	C11202 (trenutno sredstvo)
Post-Marketing Surveillance (PMS) Study FemoSeal® VCS #11202 (R1756-03) (2009)	C11202 (trenutno sredstvo)
The FemoSeal Vascular Closure System (VCS) Registry: A Prospective, Multi-Center, Observational Study in Europe (2022)	C11202 (trenutno sredstvo)
Planirane/tekuće kliničke studije nakon plasiranja na tržište	
FemoSeal™ Vascular Closure System PMCF Clinical Survey Protocol: Addressing Gaps Identified in CER-FS-2023	C11202 (trenutno sredstvo)
Aktivnosti kliničkog praćenja nakon plasiranja na tržište	
20210292 - FemoSeal Vascular Closure System Clinical Survey Report (2020)	C11202 (trenutno sredstvo)
20210298 – Guidewire Clinical Survey Report (2020)	C11202 (trenutno sredstvo)
FemoSeal™ Vascular Closure System PMCF Clinical Survey Report: Addressing Gaps Identified in CER-FS-2021 and CER-FS-2022 (2024)	C11202 (trenutno sredstvo)

Tabela 2.4 u nastavku opisuje sve studije za sredstvo FemoSeal™ VCS pre plasiranja na tržište.

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Funktionell djurtestning före marknadsintroduktion	Enhetsversion
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vasculature System. Functional Implantation Test in Sheep 12-, 15-, and 18-month Time Periods. (R1051-1), 2000, Radi Medical System AB (performed by Biomatech, France).	11200 (originalenhet)
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep 12-, 24-, and 36-Week Time-Periods. (R1050-1), 2002, Radi Medical System AB (performed by Biomatech, France).	11200 (originalenhet)
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep and in the Pig for 5 Weeks, (R1052-01) 2002, Radi Medical System AB (performed by Biomatech, France).	11200 (originalenhet)
Första kliniska studien på människa	
Multi-Centre Clinical Trial of the FemoSeal Vascular Closure System for Sealing Femoral Arterial Punctures After Diagnostic/Interventional Cardiology Procedures (EU-SE-001); 2004, Radi Medical Systems AB.	11200 (originalenhet)
Slutförda kliniska studier efter marknadsintroduktion	
PMS FemoSeal® VCS #11202 (R1756-01) (2006)	C11202 (aktuell enhet)
Confirmatory Assessment of FemoSeal™ Vascular Closure System in 7F-Sheath. (2008)	C11202 (aktuell enhet)
Post-Marketing Surveillance (PMS) Study FemoSeal® VCS #11202 (R1756-03) (2009)	C11202 (aktuell enhet)
The FemoSeal Vascular Closure System (VCS) Registry: A Prospective, Multi-Center, Observational Study in Europe (2022)	C11202 (aktuell enhet)
Planerade/pågående kliniska studier efter marknadsintroduktion	
FemoSeal™ Vascular Closure System PMCF Clinical Survey Protocol: Addressing Gaps Identified in CER-FS-2023	C11202 (aktuell enhet)
Kliniska uppföljningsaktiviteter efter marknadsintroduktion	
20210292 - FemoSeal Vascular Closure System Clinical Survey Report (2020)	C11202 (aktuell enhet)
20210298 – Guidewire Clinical Survey Report (2020)	C11202 (aktuell enhet)
FemoSeal™ Vascular Closure System PMCF Clinical Survey Report: Addressing Gaps Identified in CER-FS-2021 and CER-FS-2022 (2024)	C11202 (aktuell enhet)

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Pazarlama Öncesi Fonksiyonel Hayvan Testi	Cihaz Sürümü
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vasculature System. Functional Implantation Test in Sheep 12-, 15-, and 18-month Time Periods. (R1051-1), 2000, Radi Medical System AB (performed by Biomatech, France).	11200 (orijinal cihaz)
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep 12-, 24-, and 36-Week Time-Periods. (R1050-1), 2002, Radi Medical System AB (performed by Biomatech, France).	11200 (orijinal cihaz)
Biocompatibility and Biofunctionality of a Biodegradable Implant to Achieve Haemostasis in the Vascular System. Functional Implantation Test in the Sheep and in the Pig for 5 Weeks, (R1052-01) 2002, Radi Medical System AB (performed by Biomatech, France).	11200 (orijinal cihaz)
İnsanlarda İlk Klinik Çalışma	
Multi-Centre Clinical Trial of the FemoSeal Vascular Closure System for Sealing Femoral Arterial Punctures After Diagnostic/Interventional Cardiology Procedures (EU-SE-001); 2004, Radi Medical Systems AB.	11200 (orijinal cihaz)
Tamamlanan Pazarlama Sonrası Klinik Çalışmalar	
PMS FemoSeal® VCS #11202 (R1756-01) (2006)	C11202 (mevcut cihaz)
Confirmatory Assessment of FemoSeal™ Vascular Closure System in 7F-Sheath. (2008)	C11202 (mevcut cihaz)
Post-Marketing Surveillance (PMS) Study FemoSeal® VCS #11202 (R1756-03) (2009)	C11202 (mevcut cihaz)
The FemoSeal Vascular Closure System (VCS) Registry: A Prospective, Multi-Center, Observational Study in Europe (2022)	C11202 (mevcut cihaz)
Planlanan / Devam Eden Pazarlama Sonrası Klinik Çalışmalar	
FemoSeal™ Vascular Closure System PMCF Clinical Survey Protocol: Addressing Gaps Identified in CER-FS-2023	C11202 (mevcut cihaz)
Pazar Sonrası Klinik Takip Faaliyetleri	
20210292 - FemoSeal Vascular Closure System Clinical Survey Report (2020)	C11202 (mevcut cihaz)
20210298 – Guidewire Clinical Survey Report (2020)	C11202 (mevcut cihaz)
FemoSeal™ Vascular Closure System PMCF Clinical Survey Report: Addressing Gaps Identified in CER-FS-2021 and CER-FS-2022 (2024)	C11202 (mevcut cihaz)

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1. Alabduljabbar M, Alhamdan B, Alabdrabalnabi M. Access site complications of manual compression versus closure devices after lower limb revascularization: a systematic review. *International Journal of Medicine in Developing Countries*. 2022;6(12):1572-1580. doi:10.24911/ijmdc.51-1666716459
2. Chugh Y, Bavishi C, Mojadidi MK, et al. Safety of transradial access compared to transfemoral access with hemostatic devices (vessel plugs and suture devices) after percutaneous coronary interventions: A systematic review and meta-analysis. *Catheter Cardiovasc Interv*. Aug 2020;96(2):285-295. doi:10.1002/ccd.29061
3. Ciprian Cacuci A, Krankenberg H, Ingwersen M, et al. Access Site Complications of Peripheral Endovascular Procedures: A Large, Prospective Registry on Predictors and Consequences. *J Endovasc Ther*. Oct 2021;28(5):746-754. doi:10.1177/15266028211025044
4. Faraj J, Daneshmand A, Mwipatayi BP. Large spontaneous median sacral artery pseudoaneurysm masquerading as hip pain. Article. *J Vasc Surg Cases Innov Tech*. Jun 2020;6(2):224-227. doi:10.1016/j.jvscit.2020.02.012
5. Gabrielli R, Siani A, Smedile G, Rizzo AR, Accrocca F, Bartoli S. Carotid Artery Stenting versus Carotid Endarterectomy in Terms of Neuroprotection DW-MRI Detected and Neuropsychological Assessment Impairment. Article. *Ann Vasc Surg*. Jan 2024;98:68-74. doi:10.1016/j.avsg.2023.05.046
6. Goueffic Y, Picquet J, Schneider F, et al. A Randomized Trial Comparing Polymer Versus Suture- Based Vascular Closure Devices for Arterial Closure Following Lower-Limb Arterial Endovascular Revascularization. *Cardiovasc Intervent Radiol*. Dec 2021;44(12):1883-1892. doi:10.1007/s00270-021-02940-z
7. Ha TH, Lee SH, Park SJ, Chung HH. Could real-time ultrasonography guidance be useful for the effective deployment of FemoSeal in common femoral arteriotomy? *Ultrasonography*. Jul 2021;40(3):449-454. doi:10.14366/usg.20127
8. Hassu J, Meyborg M, Gebauer K, Makowski L, Malyar N. 51. Embolization of closure device particle as a rare case of new-onset claudication - Normal diagnostic findings despite limiting symptoms. A Diagnostic Challenge. Conference Abstract. *Vasa*. Nov 2022;51(Suppl 109):1-21. doi:10.1024/0301-1526/a001037
9. Kaszczewski P, Kozubek H, Ostrowski T, et al. Rapid Formation and Hybrid Treatment of a Large Superior Mesenteric Artery Aneurysm. Article. *Am J Case Rep*. Oct 29 2023;24:e939558. doi:10.12659/AJCR.939558

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